



Operation Manual

Coatron M4





Instrumentation and reagents for human coagulation and hemostasis
Copyright © 2014, TECO GMBH
Revision 14
Firmware 1.13d
Issue Jan-2014
Document No: 38 427 01

Warning	This manual is valid for firmware M4 1.13d. The manual may differ slightly from the actual product as a result of product improvements. Please read the Operation Manual in its entirety prior to operate. In order to ensure a high level of performance, all safety information must be followed.
Copyright	Copyright © 1999 - 2014 by TECO GmbH Neither the Operator's Manual nor any part thereof may be copied, digitally processed or otherwise transferred without written permission from TECO GmbH
Trademarks	Coatron M4 is a trademark of TECO GmbH, Germany. Other product names used in this Operator's Manual are trademarks of the respective companies.
Software	The software for TECO GmbH products is the intellectual property of TECO GmbH, which company retains all rights to usage of the software. The purchaser of a Coatron M4 acquires rights of use for this software.
Warranty	TECO GmbH guarantees that the instrument will delivered in a fault free condition. Any damages resulted by accidents , improper use, using not recommended material or negligence of maintenance are excluded from warranty. Warranty will expire if not authorized persons perform any service on the instrument
Service	Repairs to the instrument may only be carried out by trained personnel, and replacement parts must comply with instrument specifications. Please contact the local distributor of TECO GmbH if service is required.

Manufacturer	TECO GmbH Dieselstr. 1 84088 Neufahrn i.NB Germany
	Phone: +49 (0) 8773 70780 00
	Fax: +49 (0) 8773 70780 29
	Email: info@teco-gmbh.com
	Internet: http://www.teco-gmbh.com

The instrument is "Made in Germany".

Firmware History Coatron M4

- 1.01 - First Release
- 1.02 - Keyboard Check introduced in service menu
- Duplicate determination introduced
- 1.03 - Barcode functionality introduced
- 1.04 - 6 dialogue languages introduced
- Autosensitivity for clotting method
- New Fib algorithm
- 1.05 - 6 dialogue languages introduced
- Autosensitivity for clotting method
- New Fib algorithm
- Result Flag R & S introduced
- Autostart introduced
- 1.06 - Chromogenic method introduced (requires UV-LED optic upgrade)
- New tests D-Dimer , Eccarin ECAH & ECAT
- Coag & Time correction introduced
- 1.07 - Quick change with direct test keys
- fixed profile measurement PT , PTT , TT & TT
- 1.08 - Tecam Smart interface implemented (LIS connection)
- Improved optic check for hi lipemic, bilirubin or other disturbance
- Autostart sensitivity can be adjusted for every test individually
- General improvements of measurements
- D-Dimer can be adjust with time & od correction for any reagent
- 1.09 - Profile can be set free
- Bugfix Tecam Interface & Print-Out
- D-Dimer default adjustment for Amax Auto D-Dimer (400nm)
- 1.10 - Print routine optimized for new 32 character printer
- Quick test change with direct keys or up/down key
- Results (s,%,INR,...) are scrolled
- Result Flag 'B' (Biphasic aPTT) introduced
- Result Flag 'F' (low fibrinogen) introduced
- new test LA (DRVVT) with automatic ratio calculation
- new test APCR with automatic ratio calculation
- statistic counter for PT,aPTT,FIB,AT,DD and all introduced
- bugfix for high bilirubin samples
- Autostart sensitivity can be set now within SETUP TEST menu
- 1.11 - Optic check before menu analyse removed
- D_Dimer adapted to new TECO Blue D-dimer reagent 2008
- Bugfix DD: For zero D-Dimer curves the former system could reboot sometimes
- 1.12b - All factors are reported in IU/mL to display, printout and host
- Stopwatch is only visible if in use
- Optic amplification is now doubled for extreme turbid samples
- A press on key "UNIT" after measurement will:
 - x stop unit autoscrolling in display
 - x restore old result
- Maximum reading time depending method. Maxtime PT based=300s,Maxtime aPTT based=450s.
- Sensitivity of clotting algorithm slightly increased - especially at begin of reaction.
- A press to key MENU during analysis will return to main menu, if no measurement is running.
With former firmware even just activated channels had blocked the menu.
- Kinetic algorithm can now be adjusted to required time scale by hidden value "COAG Correction".

0% (Default)	means:	Endtime=100s	Starttime=Endtime/10 = 10s
100%	means:	Endtime=200s	Starttime=Endtime/10 = 20s
-20%	means:	Endtime=80s	Starttime=Endtime/10 = 8s

- Limits for units
new: %=1-199 U=1-999 (DD=1-5000) I=0.3-15 R=0.3-15
old: %=3-199 U=0-750 (DD=0-5000) I=0.5-12 R=0.5-9

-
- 1.13b - Important bugfix:
A print of mean result on channel 1+2 can cause prolonged results on consequent channels 3+4. This is fixed now.
- Barcode scanner improvement:
x Depending on barcode scanner, sometimes the label was read falsely.
x Coatron M4 can now be operated with TECAM(LIS)+Scanner at same time.
Before it was either barcode or LIS.
x Barcode label is limited to max. 9 characters
x Barcode can now be input either by scanner or manually at same time.
x A barcode reading event will now activate the next optic channel.
- New method "MECHANIC":
This method is replacing the older method "100mOD".
In principle both method are identical, only the stop limit was changed from 100mOD to 25mOD. Clotting tests (PT,PTT,...) run with method "MECHANIC" (instead of CLOT) will give shorter second results and can be better correlate to mechanical instrument
- TECAM/LIS improvement:
Zero results (eg. D-Dimer= 0ng/mL) were treated before always as no result. Now the Coatron M4 will send differently. As zero without any units, if the result is invalid, otherwise as "<1" with units.
(Example valid zero for DD: 0mE 0ng/mL-> transmission= 0mE <1ng/mL)
(Example invalid zero for PT: sec=+++ 0%-> transmission= 0 0)
- Measurement improvement:
x Better mixing: With old version, the optic was calibrated immediately after start. During this procedure it is not allowed to touch or mix the cuvette. For D-Dimer testing mixing is explicitly required.
Now the system calibrate short before end of deadtime. This gives operator much better time window for mixing.
Additionally a countdown gives information about time to mix.
- x Life time of LED
The optic LEDs are only powered during measurement.
This will increase life time of LEDs.
- x Calibration mode with key "5"
All channels can now be activated in duplicate mode with key "5".
The difference to key "0" or key "." activation is, that the optic is switched to high sensitive mode. This will allow to measure also highest dilutions, like used during calibration.
High sensitive mode is displayed with blinking "CALIBRATE".
- Scroll of units synchronized
Instrument displays the results (s,%,INR,..) in an blinking mode.
Sometimes this can be confusing. A press on any of the UNIT keys will stop blinking for all channels and synchronize the units.
(eg. report all channels in seconds)
- 1.13c - Important bugfix:
"+++" results caused a system reboot. This is fixed now.
- 1.13d - TIMER ALARM feature:
It can be set to OFF, 30sec, 60sec, 500sec.
OFF = TIMER keys causes regular stopwatch.
>0s = TIMER keys causes a countdown
- Important bugfix:
Profile testing (STAT): Results were falsely transmitted to TECAM

TABLE OF CONTENTS

1.0	SAFETY INFORMATION	10
2.0	GENERAL	11
2.1.	Intended purpose	12
2.2.	Installation	12
2.3.	Technical data	13
2.4.	CE marking	14
3.0	INSTRUMENT COMPONENTS	15
3.1.	Incubator Block	15
3.2.	Control panel	16
3.3.	Rear of equipment	17
4.0	THEORY OF OPERATION	18
4.1.	Clotting Assay (CLOT)	19
4.2.	Mechanical clotting time (MECHANIC)	19
4.3.	Derived Fibrinogen (CLOT + FIB)	19
4.4.	Chromogenic Assay (KINETIC)	20
4.5.	Chromogenic Ecarin Time (MECHANIC)	20
4.6.	Immunturbidimetric Assay (IMMUNO)	21
5.0	OPERATING INSTRUCTIONS	22
5.1.	Setup System	22
5.1.1.	Language	22
5.1.2.	Fibrinogen Concentration Units	22
5.1.3.	Temperature Control	23
5.1.4.	Signal	23
5.1.5.	Autostart	23
5.1.6.	Timer Alarmt	23
5.1.7.	Contrast of the LCD (Liquid Crystal Display)	23
5.1.8.	Speed of the Mixer	23
5.1.9.	Patient Identification	24
5.1.10.	Date/Time	24
5.2.	Setup Test	26
5.2.1.	Setup Test	26
5.2.2.	Units	27
5.2.3.	Standard Curve	27
5.2.4.	Correlation Factor (linearity index for calibration data)	28
5.2.5.	Store Data	28
5.2.6.	Print Test	28
5.2.7.	Autostart	28

5.3. Test Analysis	30
5.3.1. Test Selection	30
5.3.2. Optic Activation	31
5.3.3. Entering Patient Identification Numbers (PID)	31
5.3.4. Duplicate testing	32
5.3.5. Starting the Analysis	32
5.3.6. Display during measuring	33
5.3.7. Manual break of measurement	33
5.3.8. Return to main menu	33
5.3.9. Unit Key Functions	33
5.3.10. Stopwatch Functions	33
5.3.11. Alarm Functions	33
5.3.12. Result Warning Messages	34
5.4. Hidden Functions	35
5.4.1. Set system to default	35
5.4.2. Change temperature adjustment	36
5.4.3. Set all test calibration points to zero	36
5.4.4. Set OD correction	36
5.4.5. Set COAG correction	37
6.0 SERVICE MENU	38
6.1. System Analysis	38
6.2. Optic-Values	39
7.0 TROUBLESHOOTING	41
8.0 MAINTENANCE	44
8.1. Recommended Maintenance	44
8.2. Temperature Adjustment	44
8.3. Cleaning procedures	44
9.0 APPLICATIONS	45
9.1. Prothrombin Time	46
9.2. Derived Fibrinogen	47
9.3. Clauss Fibrinogen Assay	48
9.4. Thrombin Time Assay	49
9.5. APTT	50
9.6. PT-Based Factor Assays (II, V, VII & X)	51
9.7. APTT-Based Factor Assays (VIII, IX, XI & XII)	52
10.0 SPECIAL FUNCTIONS	53
10.1. Software Upgrading	53

11.0	TECAM SMART THE LIMS SOLUTION	54
11.1.	General	54
11.2.	Interface Protocol	55
12.0	PRODUCT CATALOGUE	58

TABLE OF FIGURES

Figure 1 - Incubator Block	15
Figure 2 - Control Panel	16
Figure 3 - Rear of Equipment	17
Figure 4 - The detection principle	18
Figure 5 - The Turbidity Method	19
Figure 6 - The chromogenic method	20
Figure 7 - Latex agglutination	21
Figure 8 - Relationship of light absorbance and concentration of D-dimer	21
Figure 9 - Flow diagram for the "Setup" System" Submenu	25
Figure 10 - Flow Diagram for the "Setup Test" Submenu	29
Figure 12 - Flow Diagram of "SERVICE" Submenu	40
Figure 13 – LIMS communication	55

Symbols

Symbol	Meaning	Explanation
	Advice	Indicates important information and tips.
	Warning!	Risk of possible health damage or considerable damage to equipment if warning is not heeded.
	Biohazard !	Equipment can be potentially infectious due to the samples and reagents used .
	Danger!	Potential risk to operating personnel or equipment due to electric shock.

1.0 Safety information



Recommend materials

Use only origin disposables.
Use only manufacturer approved material.



Avoid contact

Never touch moving parts such as the measurement rotor or pipetting arm during device operation.



Do quality control

Carry out control measurement runs at regular intervals to ensure that the Analyzer continues to function faultlessly.



Waste cuvettes

The cuvette blocks are intended as single-use items only.



Infectious Material

Avoid direct contact with samples and sample residues in the used cuvettes.

Infectious material such as cuvette waste and liquid waste must be disposed in compliance with local regulations governing for infectious materials.

Wear medical infection grade protective gloves for all cleaning and maintenance work involving potential contact with infectious liquids and use each pair of gloves once only.

Use a hand disinfectant product, e.g. Sterilium®, to disinfect your hands after completion of the work.



Environmental condition

Ambient temperature must be 18 – 25°C
humidity must be below 80%
avoid any vibrations or impacts to analyzer
do not use analyser if explosive or inflammable gas is around.



Elektrical Safety

Make sure the operating voltage setting is correct before connecting the device to the power mains.
Use only shockproof (grounded) electrical sockets.
Use only shockproof extension leads in perfect condition. Defective leads must be replaced without delay.
Never intentionally interrupt protective ground contacts.
Never remove housing elements, protective covers or secured structural elements, since so doing could expose parts carrying electric current.
Make sure surfaces such as the floor and workbench are not moist while work is being done on the device.



Safety Information for operation

If the instrument is used in a manner not specified by the manufacturer, the protection impairment could be affected!

2.0 General

The **Coatron M4** is a manual 4 channel photo-optical instrument that offers clotting & chromogenic & immunturbidimetric testing capabilities. The **Coatron M4** can be used for a wide variety of coagulation and fibrinolysis tests such as:

• Prothrombin Time (PT) • Activated Prothrombin Time (aPTT) • Thrombin Time (TT) • Venom Time (VT) • Fibrinogen (FIB) • Factors (FII - FXII) • Antithrombin (AT3) • Heparin (HEP) • Activated PC resistance (APCR) • Lupus Anticoagulant (Screen, Confirm)	• Protein-C (PC) • Protein-S (PS) • D-Dimer (DD) • von Willebrand Factor (VWF) • Ecarin Chromogenic Assay Thrombin (ECAT) • Ecarin Chromogenic Assay Hirudin (ECAH) • Plasminogen (PLG) • a2-Antiplasmin (A2AP)
---	--

FEATURES:

- Coagulation analyser for turbidimetric, chromogenic and immunturbidimetric assays.
- Highly reliable, longlifying and nearly service free system
- Autosense optics to eliminate interferences like Bilirubin, Hemoglobin
- Approved clotting algorithm for all kind of samples and reagents. If there is a clot - it will be detected.
- Biphasic aPTT waveform curve detection to DIC indication
- Low fibrinogen curve detection
- Fibrinogen concentration can be derived from a PT result additionally. For sure the standard CLAUS method is also available.
- Calculation in Activity %, INR, Ratio, g/L or mg/dL
- Every test is programmable up to 5 calibration points
- Correlation analyse of calibration curve
- 4 Stop-watch functions which can be used independently
- 4 Alarm timer functions which can be used independently
- Multi-language software
(German, English, Spanish, French, Italian, Portuguese)
- Patient identification (No PID, manual input, autoseries, barcode)
- Test-Duplicate Modus
- Free Profile testing (PT, APTT, FIB)
- APC-R with automatic ratio calculation
- DRVVT with automatic ratio calculation
- Micro volume testing (total of 60 - 75µL)
- Reagent stirring with magnetic bars
- Routines for selftest (trouble-shooting)
- Routines for print outs (result, calibration, service, system)
- Optional autopipette for electronic triggered start
- Automatic Start triggered by adding reagent
- Optional data management and research software
- Optional printer
- Optional barcode scanner for patient identification
- Easy software update
- Interface for Laboratory Information & Management Systems (LIMS)
- CE marked
- Small dimension and weight - fits on every desktop

2.1. Intended purpose

The **Coatron M4** is designed to carry out coagulometric tests such as PT, PTT, TT, fibrinogen, single factor tests, chromogenic tests, like Antithrombin-III and immunturbidimetric tests, such as D-dimer).

Use only citrate plasma for test analysis runs: Mix 9 parts venous blood with 1 part 3.2% (0.105M) sodium citrate and centrifuge the mixture at 1500g for approx. 10 minutes. Plasma must be used within 4h.



Do not use plasma with more than 25mg/dL Bilirubin concentration
Do not use plasma with more than 1000mg/L Hemoglobin concentration

The **Coatron M4** must be operated by a specialist trained in clinical laboratory techniques who has also received instruction and training in operation of the **Coatron M4** and has read and understood this Operator's Manual

2.2. Installation

No special precautions are necessary when starting up the **Coatron M4**. However, the following is recommended:

- Place on a level surface in an area free from excessive temperature fluctuations
- Avoid vibration during measurement
- Protect the instrument from direct sunlight, moisture and dust
- Check that the voltage and frequency data on the identification plate on the instrument agree with the local power rating before starting the instrument for the first time



The instrument is connected to the power supply by the main cable supplied.
If obvious damage has occurred during shipping, **do not use**.
Contact your local distributor for replacement or repair.



The power cord always has to be easily accessible during normal operation.

2.3. Technical data

Instrument:

Boards	SMD (Small Mounted Devices) based
Microprocessor	NEC V25
Flash-EEPROM	128 KByte
RAM	128 KByte
EEPROM	2 KByte
AD-Converter	18 Bit (16 bit used)
Optics	4 LED's ultra bright, pulse modulation control
RS 232	9600 Bauds, 8 Data, 1 Stop, no parity

Dimensions:	LxWxH: 280 x 225 x 115 mm
Weight (net)	2,3 kg

Performance:	Precision: CV <5%
--------------	-------------------

Power Supply:

external, 42 W max	
Input voltages	100 Vac to 240 Vac / 47 to 63 Hz
Output voltages	+5Vdc/5A; +15Vdc/2A; -15Vdc/0,8A
Approvals	CE, EN60601-1, UL, RoHS

Keyboard:

4x8 matrix, foil keyboard, with Test, Function and numerical keys

Display:

4 lines x 20 characters Liquid Crystal Display

Incubation block:

6x4 double-cuvette prewarming positions, 4 measuring and 5 reagent positions

Ambient conditions:

Operating Temperature	+15 °C to +30 °C
Humidity	< 70% rel. humidity
Amplitude (sea level)	< 3,000m
Free of dust	Grade 2
Impact resistance	according to IEC/EN 61010-1, 8.2.2
Not allowed	Vibrations, direct sun light; direct exposure of air condition.

Storage and Transport	Temperature range: 0 °C to +40 °C
-----------------------	-----------------------------------

2.4. CE marking (example)

	TECO KONFORMITÄTSERKLÄRUNG DECLARATION OF CONFORMITY Wir / We TECO Medical Instruments Production and Trading GmbH Name des Anbieters / Supplier's name Dieselstrasse 1, D-84088 Neufahrn NB Anschrift / Address erklären in alleiniger Verantwortung, dass das Produkt declare under our own responsibility, that the product Coatron M4 Bezeichnung, Typ oder Modellname / name, type or model <table border="0" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; vertical-align: top;"> den Anforderungen der folgenden Richtlinien entspricht: 1. Richtlinie 98/79/EG über In-vitro Diagnostika 2. Richtlinie 2004/108/EG über Elektromagnetische Verträglichkeit </td> <td style="width: 50%; vertical-align: top;"> corresponds to the requirements of: 1. Directive 98/79/EC on In-vitro diagnostic medical devices 2. Directive 2004/108/EC on electromagnetic Compatibility </td> </tr> </table> <table border="0" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; vertical-align: top;"> Zur Beurteilung der Konformität wurden folgende harmonisierte Normen herangezogen: 1. Sicherheit: EN 61010-1:2002 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Allgemeine Anforderungen EN 61010-2-101:2003 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Besondere Anforderungen an In-vitro-Diagnostik (IVD) Medizingeräte 2. EMV: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Elektromagnetische Verträglichkeit – Anforderungen 3. Risikomanagement: DIN EN ISO 14971:2001: Medizinprodukte – Anwendung des Risikomanagement auf Medizinprodukte 4. Informationen: DIN EN 1041:1998: Bereitstellung von Informationen durch den Hersteller eines Medizinproduktes </td> <td style="width: 50%; vertical-align: top;"> The following harmonized standards have been used: 1. Safety: EN 61010-1:2002 Safety requirements for electrical equipment for measurement, control and laboratory use: General requirements for safety EN 61010-2-101:2003 Safety requirements for electrical equipment for measurement, control and laboratory use: Particular requirements for in-vitro diagnostic (IVD) 2. EMC: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Electromagnetic compatibility – Requirements 3. Risk management: DIN EN ISO 14971:2001: Medical devices – Application of risk management to medical devices 4. Information: DIN EN 1041:1998: Information supplied by the manufacturer with medical devices </td> </tr> </table> <table border="0" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; vertical-align: top;"> Das QM-System des Herstellers ist zertifiziert durch DEKRA-ITS, Stuttgart, nach: EN ISO 9001:2000 EN ISO 13485:2003/Richtlinie 98/79/EG </td> <td style="width: 50%; vertical-align: top;"> The QM-system of the manufacturer is certified from DEKRA-ITS, Stuttgart, for: EN ISO 9001:2000 EN ISO 13485:2003/Directive 98/79/EG </td> </tr> </table> Neufahrn, 05.11.2007 Ort und Datum Place and date of issue  Christian Hötzel Geschäftsführer Christian Hötzel, General Manager	den Anforderungen der folgenden Richtlinien entspricht: 1. Richtlinie 98/79/EG über In-vitro Diagnostika 2. Richtlinie 2004/108/EG über Elektromagnetische Verträglichkeit	corresponds to the requirements of: 1. Directive 98/79/EC on In-vitro diagnostic medical devices 2. Directive 2004/108/EC on electromagnetic Compatibility	Zur Beurteilung der Konformität wurden folgende harmonisierte Normen herangezogen: 1. Sicherheit: EN 61010-1:2002 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Allgemeine Anforderungen EN 61010-2-101:2003 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Besondere Anforderungen an In-vitro-Diagnostik (IVD) Medizingeräte 2. EMV: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Elektromagnetische Verträglichkeit – Anforderungen 3. Risikomanagement: DIN EN ISO 14971:2001: Medizinprodukte – Anwendung des Risikomanagement auf Medizinprodukte 4. Informationen: DIN EN 1041:1998: Bereitstellung von Informationen durch den Hersteller eines Medizinproduktes	The following harmonized standards have been used: 1. Safety: EN 61010-1:2002 Safety requirements for electrical equipment for measurement, control and laboratory use: General requirements for safety EN 61010-2-101:2003 Safety requirements for electrical equipment for measurement, control and laboratory use: Particular requirements for in-vitro diagnostic (IVD) 2. EMC: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Electromagnetic compatibility – Requirements 3. Risk management: DIN EN ISO 14971:2001: Medical devices – Application of risk management to medical devices 4. Information: DIN EN 1041:1998: Information supplied by the manufacturer with medical devices	Das QM-System des Herstellers ist zertifiziert durch DEKRA-ITS, Stuttgart, nach: EN ISO 9001:2000 EN ISO 13485:2003/Richtlinie 98/79/EG	The QM-system of the manufacturer is certified from DEKRA-ITS, Stuttgart, for: EN ISO 9001:2000 EN ISO 13485:2003/Directive 98/79/EG	
den Anforderungen der folgenden Richtlinien entspricht: 1. Richtlinie 98/79/EG über In-vitro Diagnostika 2. Richtlinie 2004/108/EG über Elektromagnetische Verträglichkeit	corresponds to the requirements of: 1. Directive 98/79/EC on In-vitro diagnostic medical devices 2. Directive 2004/108/EC on electromagnetic Compatibility							
Zur Beurteilung der Konformität wurden folgende harmonisierte Normen herangezogen: 1. Sicherheit: EN 61010-1:2002 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Allgemeine Anforderungen EN 61010-2-101:2003 Sicherheitsbestimmungen für elektrische Mess-, Steuer-, Regel- und Laborgeräte: Besondere Anforderungen an In-vitro-Diagnostik (IVD) Medizingeräte 2. EMV: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Elektromagnetische Verträglichkeit – Anforderungen 3. Risikomanagement: DIN EN ISO 14971:2001: Medizinprodukte – Anwendung des Risikomanagement auf Medizinprodukte 4. Informationen: DIN EN 1041:1998: Bereitstellung von Informationen durch den Hersteller eines Medizinproduktes	The following harmonized standards have been used: 1. Safety: EN 61010-1:2002 Safety requirements for electrical equipment for measurement, control and laboratory use: General requirements for safety EN 61010-2-101:2003 Safety requirements for electrical equipment for measurement, control and laboratory use: Particular requirements for in-vitro diagnostic (IVD) 2. EMC: EN 61326-1:1997+A1:1998+A2:2001+A3:2003 Electromagnetic compatibility – Requirements 3. Risk management: DIN EN ISO 14971:2001: Medical devices – Application of risk management to medical devices 4. Information: DIN EN 1041:1998: Information supplied by the manufacturer with medical devices							
Das QM-System des Herstellers ist zertifiziert durch DEKRA-ITS, Stuttgart, nach: EN ISO 9001:2000 EN ISO 13485:2003/Richtlinie 98/79/EG	The QM-system of the manufacturer is certified from DEKRA-ITS, Stuttgart, for: EN ISO 9001:2000 EN ISO 13485:2003/Directive 98/79/EG							

3.0 Instrument Components

3.1. Incubator Block

The incubator block is made from aluminum, which ensures equal distribution of heat. The temperature of the incubator block is regulated to 36.5°C - 37.5°C





Figure 1 - Incubator Block

3.2. Control panel



Temp.	indicates Temperature is in range (36-38°C)
Optic 1 - 4	activate/start/ stop channel 1-4
Unit 1 - 4	convert result into other dimensions
Timer 1 - 4	activate/reset stopwatch 1-4
PT, PTT, TT, STAT	direct access to test
Menu	Go back to main menu or next entry or exit menu.
Cursor up/down	Scroll up/down during - Test selection - Parameter selection - Patient-ID increment
Numeric keys	for input of any numbers
	All optic channels can be activated with key - 5/CAL: Calibration mode - 0/1x: Single mode - . /2x: Duplicate mode
Enter	confirm entry, jump to next entry

Figure 2 - Control Panel

3.3. Rear of equipment




Figure 3 - Rear of Equipment

On-Off: Standby Switch

DC-Input: Two alternative connectors to Power Supply.

Pipette: For connector to Autopipette (Cat.No. 19 029 00). Optional accessory for automatic test start. The pipette supports four different volumes (25, 50, 100 and 200 μ L)

Computer: For connector to PC (Firmware upgrading, TECAM Smart, LIMS) or connector to barcode scanner. By using an Y-adaptor this port can be used for both LIMS and BARCODE.

Printer: For connector to Thermal Printer (Cat.No. 80 535 00)



The serial ports are always set always
9600 Baud , 8 Databits , 1 Stopbit , No parity.



Warning: The barcode-scanner is powered with 5V over PIN 9 of the RS232 Interface of the analyser. Do only use scanners with that feature.

4.0 THEORY OF OPERATION

The **Coatron M4** is a highly sensitive 4-channel-photometer. A very bright LED-Optic ensures accurate and precise results, even when iteric or lipemic samples are used. The receiver signal is detected and converted to an electrical current. During the actual test the system is searching for the best signal amplification, therefore it will support a wide range of different reagents (i.e. very turbid thromboplastins or very clear reagents). Additionally the software is based on optical density (extinction), which absorbs outside light effects.




Figure 4 - The detection principle

4.1. Clotting Assay (CLOT)

The thrombin catalyzed conversion of fibrinogen to fibrin is the final reaction in the 'coagulation cascade'. Fibrin formation results in an increase in sample turbidity which is detected by the photometer. Photometric detection is started manually by pressing the "Optic" key with simultaneous addition of the test reagent. Alternatively, the reaction is started by the addition of the test reagent using the Autopipette. The time between the start of the photometric detection, and the turning point of the reaction curve (see Figure 4) is the result. The result is displayed in seconds on the Liquid Crystal Display (and printed automatically to the optional Thermal Printer.)




Figure 5 - The Turbidity Method

The diagram is representative of a typical PT curve with normal control plasma and a curve with biphasic reaction. Biphasic aPTT reaction are highly indicators of disseminated intravascular coagulation (DIC)

4.2. Mechanical clotting time (MECHANIC)

Clotting event is not defined at turn-point of reaction, but at the begin of fibrin formation (see figure above). Therefore the results are around 20% shorter than with method "CLOT"

4.3. Derived Fibrinogen (CLOT + FIB)

The derived fibrinogen is determined using the clotting method described in section above. The fibrinogen concentration in the sample is proportional to the change in optical density in the cuvette, that accompanies the conversion of fibrinogen to fibrin at the end of the reaction. The result is expressed as "E", which represents the optical density at the end-point.

4.4. Chromogenic Assay (KINETIC)

In this method, the end result is determined from the rate of optical density change. Test plasma is pre-incubated with an enzyme (i.e. - Thrombin for determination of AT-III) and residual enzymatic activity is detected by the addition of a chromogenic substrate. The concentration of the analyte in the test plasma is directly or indirectly proportional (depending on the reagent system) to the rate of substrate hydrolysis, and is reported as the mean slope of optical density per minute ($\Delta OD(E)/min$).




Figure 6 - The chromogenic method

4.5. Chromogenic Ecarin Time (MECHANIC)

The measurement principle is similar to a regular chromogenic assays. But the result is the time between start of test and when signal broke through 25mOD.

4.6. Immunturbidimetric Assay (IMMUNO)

Intensive light is able to penetrate turbid solutions, such as latex suspensions used for the determination of D-dimer concentration. Latex particles, designed specifically for automated Dimex testing, are coated with a monoclonal antibody specific for D-dimer. If D-dimer antigen is present in the sample, an antigen-antibody reaction occurs, with a simultaneously change in light transmission. The concentration of D-dimer in the sample is directly proportional to the rate of the antigen-antibody reaction. The result is reported as the mean slope of optical density per minute ($\Delta OD(E)/min$, E = Extinction, a unit of light-absorbance). The following diagram illustrates the measurement principle of the Dimex Jr.




Figure 7 - Latex agglutination

The D-dimer concentration is proportional to the rate of change in optical density. The instrument calculates the average slope of reaction, using the linear portion of the curve only.




Figure 8 - Relationship of light absorbance and concentration of D-dimer

The kinetic algorithm for D-dimer testing is illustrated with three typical reaction curves. At high doses the linear relationship between signal and concentration is not valid. This is called "High Dose Hook Effect".

5.0 Operating Instructions

This section provides general instructions necessary for the user to achieve maximal use and benefit from the **Coatron M4**. Please refer to section 9.0 for specific test applications.

The On/Off switch is located on the rear panel of the instrument. **For optimal results, do not operate until the temperature indicator light is on.** It takes approximately 10-15 minutes for the instrument to equilibrate to 37°C. The general sequence of operation for test analysis is:

- 1.) Enter the "SETUP SYSTEM" submenu to confirm/change system settings
- 2.) Enter "SETUP TEST" to select test parameters and enter calibration data if desired; and
- 3.) Enter the "ANALYSIS" submenu for sample testing.

From the **Main Menu**, the following options are available:

1. ANALYSIS
 2. SETUP TEST
 3. SETUP SYSTEM
 4. SERVICE

At each screen, selections are made using the Up/Down cursors. To proceed to the next menu item, press either the "Menu" or "Enter" key. If a mistake is made, press the "Menu" key until the main menu appears and start over.



**To return the system to default values,
Press simultaneously "Optic 1" + "." + "Enter" keys.**

5.1. Setup System

To enter this submenu, press #3 from the main menu. The default values for the system parameters are:

LANGUAGE:	ENGLISH
FIBRINOGEN:	mg/dL
TEMP.CONTROL:	ON
SIGNAL:	ON; VALUE 325
AUTOSTART:	ON
TIMER ALARM:	OFF
CONTRAST OF LCD:	VALUE: 25
SPEED OF MIXER:	VALUE: 215
PAT.IDENT.:	NO PAT.ID.
INTERFACE:	9600 Baud 8Data 1Stop No parity

5.1.1. Language

English, German, French, Italian, Spanish, Portuguese
Use the cursor keys to select, press "Enter" or "Menu" to advance

5.1.2. Fibrinogen Concentration Units

Use the cursor keys to select mg/dL or g/L, press "Enter" or "Menu" to advance

5.1.3. Temperature Control

On/Off - use the cursor keys to select, press "Enter" or "Menu" to advance. For temperature adjustment, refer to section **8.2 Temperature Adjustment**

5.1.4. Signal

ON/OFF (A beep at the start and end of testing) - use the cursor keys to select, press "Enter" or "Menu" to advance.

Higher/lower - use the cursor keys to change and "Enter" to advance.

5.1.5. Autostart

Change with cursor keys to change and continue with "ENTER"

ON Measurement starts with adding of reagent automatically.
No need for use of the Autopipette or extra pressing of the relevant Optic channel key.

Off regular mode for start with optional available Autopipette or manual start
by pressing relevant Optic channel key



Activate optic channel just before adding of reagent. Movements of the cuvette can pre-start measurement – Do not touch cuvette, if optic channel is active !



The sensitivity of the autostart can set for every test individually. Refer to section 5.2.7 Autostart

5.1.6. Timer Alarm

Change with cursor keys to change and continue with "ENTER"

OFF TIMER keys will start , stop stopwatch function

ON Set to 30s upto 500s. TIMER keys will start , stop countdown function with alarm

5.1.7. Contrast of the LCD (Liquid Crystal Display)

Higher/lower - use the cursor keys to change and "Enter" to advance.

5.1.8. Speed of the Mixer

Higher/lower - use the cursor keys to change and "Enter" to advance.

5.1.9. Patient Identification

Four choices are available:

- No Patient ID
- Autoseries
- Extra Input / Barcode
- Autoseries

Use the cursor keys to select, press "Enter" or "Menu" to advance.

No Patient ID: Results will be printed out without a patient identification number.

Autoseries Sample-id is automatically incremented by one.
Series will start at entered value (e.g. 1000, 1001, 1002, ...)

Extra Input/Barcode: Patient-Id can be entered manually or by barcode (max. 9characters)

5.1.10. Date/Time

U.S. date format (mm/dd/yy) will be selected if "English" is chosen as the language. The date format is dd/mm/yy if German is selected. The Time format is not changeable.
To return to the main menu, press "Menu".




Figure 9 - Flow diagram for the "Setup" System" Submenu

5.2. Setup Test

The specific parameters for each test are entered in this menu. Once the initial data is entered and saved it is not necessary to go into this menu for routine testing. The submenu of "**SETUP TEST**" is illustrated in Figure 6. To enter this submenu, press #2 from the main menu. The default values for "**SETUP TEST**" are:

METHOD:	CLOTTING (or Clauss for FIB)
UNITS:	seconds
CALIBRATION CURVE:	Reset to zero
AUTOSTART:	500 (250 for FIB)

5.2.1. Setup Test

To select a test by enter the numeric code of the designated test. Alternatively, the Up/Down arrow keys can be used to scroll through the entire test menu. For example, key in #1 to select PT. Press "Enter" to confirm selection. The "METHOD" can also changed when blinking.

TEST:	PT
METHOD:	CLOT
AUTOSTART:	800

Depending which test is active, following methods can be selected:

- | | |
|------------|---|
| • CLOT | Clotting assays |
| • CLOT+FIB | Fibrinogen will be derived from PT |
| • CLAUSS | Fibrinogen according to CLAUSS |
| • KINETIC | Chromogenic Assay |
| • MECHANIC | Mechanical clotting time or Chromogenic Ecarin time |
| • IMMUNO | Immunturbidimetric Assays |

5.2.2. Units

Every result is displayed in seconds [s]. However, the user can also choose to display PT results in % (% activity), R (ratio) and I (INR). Calibration data, the mean normal PT value, and/or the ISI of the thromboplastin reagent must be entered to obtain results in %, R and I. Refer to next section for information on calibration data entry. Use the Up/Down cursor keys to select the desired units, "Enter" to confirm and "Menu" to advance.

UNITS:	PT
(s-%-R-I):	s-%-I- -
NORMALVALUE:	12.2s
ISI-VALUE:	1,05

The **Coatron M4** reports results using the following units (which are test dependent):

E = Extinction (optical density)	precision X.XXX
s = seconds	precision XXX.X
R = ratio	precision XX.XX
I = INR	precision XX.XX
% = percent activity	precision XXX.XX
U = mg/L (except FIB: mg/dL or g/L)	precision XXX.X

5.2.3. Standard Curve

To obtain results in units of concentration (mg/dL, IU/mL...) or % activity, a calibration curve is needed. A minimum of **two** points is required, with a maximum of five. Calibration data is obtained by testing plasma [in duplicate (2) or quadruplicate (4)] in the "Analysis" mode. An example of calibration data entry is shown below.

Example: A PT calibration curve with derived fibrinogen. Two different calibration curves are required. The order of entry is not critical, the instrument will automatically sort calibration points.

3-point Calibration Curve for PT

100% = 12.2s
50% = 18.0s
25% = 27.2s
0% = 0.0s
0% = 0.0s

4-point Calibration Curve for Derived Fib.

591 mg/dL = 0.413E
377 mg/dL = 0.246E
267 mg/dL = 0.140E
95 mg/dL = 0.042E
0 mg/dL = 0.0E



A correct calibration curve is required to obtain results in units of concentration or % activity. Those points remaining with no data entry are not used in the calibration curve calculation. For those tests that require a zero calibration point a value greater than zero must be entered (i.e. 0,1 % and 0.1s). All calibration points can be reset to zero by simultaneously pressing the "0" and "Enter" keys. A full parameter reset will eliminate calibration data for all assays.

5.2.4. Correlation Factor (linearity index for calibration data)

Insert correct calibration data and confirm storing and printing. All test-parameters and in addition a correlation factor will be printed. The correlation factor (R^2) indicates the linearity of the calibration curve. It is 1.000 if the points are exactly on one line.

If R^2 is less than 0.980, a calibration curve with more than 3 calibration points is recommend.

5.2.5. Store Data

Yes/No. Use the cursor key to select.

5.2.6. Print Test

Yes/No. Use the cursor key to select.

5.2.7. Autostart

The sensitivity of the autostart feature can be adjusted for every test individually.

The value represents the required optical change before the instrument triggers the measurement start.

Range of sensitivity:

- Very sensitive: 300 – 500
- Normal: 500 – 1000
- Insensitive >1000
- The default value is 500.
- Increase value if test does start before adding the reagent
- Decrease value if test does not start at all.
- Set to 0 will disable the autostart for the test

Some tips:


- Pipeting clear in clear suspension will produce only low optical change and require a sensitive setup. (example: FIB or chromogenic tests)
- Try to pipet direct into the suspension instead of cuvette walls.
- Try to pipet with high pressure



Do not use autostart below 250 !

The optic could be triggered just by touching the instrument.

Figure 10 - Flow Diagram for the "Setup Test" Submenu




5.3. Test Analysis

To enter the submenu "Analysis", press #1 from the main menu and select required test with cursor keys or direct test code (e.g. PT=01) or direct keys (PT/PTT/TT/STAT). The direct test code are labeled on instrument.

5.3.1. Test Selection

The first screen in the "Analysis" submenu will ask the user to make a test selection. This can be done three ways: by pressing one of the 4 direct test keys; scrolling through the test menu with the cursor Up/Down keys; or by entering a numeric test code (i.e. - "01" for PT test). The test selection is confirmed with the **"Enter"** key.



STAT / Profile:

All 4 channels can be set individually, If test "STAT" is selected. Press the corresponding optic key and change the test with cursors keys. Confirm test by pressing optic key again.



Within the analysis menu the actual test can switched by pressing the direct test keys or cursors keys if no measurement is ongoing.

5.3.2. Optic Activation

PT: <i>ACTIVE</i>	01:30
PT: ACTIVE	00:00
PT: ACTIVE	00:00
PT: ACTIVE	00:00

The blinking message "ACTIVE" will indicate a ready channel.

The optic can be activated in four different ways.

Single activation :

- Press optic key
 - Enter / Change PID
 - Press optic key again → PID is confirmed and channel active
- Repeat procedure for each desired optic channel

Multi activation with single PID

- Press key "." (dot)
 - All four channels will activated
- All patient ID number will set logical if autoseries is active (PID = 100,101,102,103)

Multi activation with double PID

- Press key "0"
- All four channels will activated in duplicate mode

Multi activation for test calibration

- Press key "5"
- All four channels will activated in duplicate mode with highest clotting sensitivity. This is especially required for running very high diluted samples like during calibration.

5.3.3. Entering Patient Identification Numbers (PID)

The entry of a PID will be always entered in following ways

- Press key "OPTIC" to select channel
- Enter or change PID
- Press key "OPTIC" again to confirm PID. Optic is now active and ready for start.

The different patient-id entry mode can be set within menu SETUP SYSTEM.

"No Pat Id" will disable PID. Result will be printed without any PID. Press Optic 1 to confirm and activate channel. Repeat for remaining channels.

"Autoseries" will be automatically incremented PID by one for each new measurement. Press the "Optic" key and PID will be displayed. Change number optional with cursor keys. Press "Optic" again to confirm PID and repeat procedure for remaining channels.

"Manual Input/Barcode " enables individual input of PID. Press the "Optic" key and last entered number will be shown. Change PID with cursor keys or overwrite with numeric keys or read barcode. Press "Optic" again to confirm PID and repeat procedure for remaining channels. A barcode scanner event will automatically activate the next channel.

5.3.4. Duplicate testing

If duplicate testing is desired, enter the same patient identification number for channels 1 & 2 or 3 & 4 twice. The mean result will automatically be printed along with the individual channel data.



If the duplicates differ more than 7.5% from the mean value, the mean result will be flagged with 'X' !



Pressing key "0" during optic activation, will activate all channels in duplicate mode

5.3.5. Starting the Analysis

The analysis can be started, if optic is "ACTIVE" (see chapter above) by one of following events

- Press on key "OPTIC"
- Optical Autostart (pipeting reagent with normal pipet)
- Autopipette (using a pipet with an electronic trigger signal)

Step by step procedure for PT measurement

Place the required number of cuvettes in the cuvette pre-warming positions and pipet 50µL plasma in each cuvette. Press the "Timer 1" button to start the stopwatch and wait 1-10min before continue with next step. Transfer cuvette to the measurement position. Press key "OPTIC" until channel is active.

PT: <i>ACTIVE</i>	01:30
PT: ACTIVE	00:00
PT: ACTIVE	00:00
PT: ACTIVE	00:00

Now 100 uL of prewarmed Thromboplastin reagent into the first left cuvette position. If autostart is enabled, the optic will start automatically. A beeping noise will indicate the start of the reaction. Repeat for the remaining cuvettes. The reaction can also be started either by pressing the optic key or by using the Autopipette, which electronically triggers the optic channel.

In summary, press the optic key once for patient identification data entry, press again to activate the channel, and a third time to manually initiate the reaction. (Do not press a third time if using the Autopipette to start the reaction!)



If using the Autopipette,

- **pipet always from left to right (channel 1 - 4)!**
- **disable the autostart function in the SETUP SYSTEM**

5.3.6. Display during measuring

Once started, a short beeping noise is followed by a countdown. During the countdown it is allowed to mix or touch the cuvette. After the countdown the actual optical density will be displayed. Avoid contact with the cuvette while this message is shown. A beeping noise will sound when the reaction is complete and the result will be displayed on the screen. If the thermal printer is attached, results will automatically be printed.

PT:	050	00:00
PT:		00:00
PT:		00:00
PT:		01:00

Countdown after test start.
Touching or mixing cuvette
is allowed.

PT:	000	00:00
PT:		00:00
PT:		00:00
PT:		01:00

Countdown finished.
Don't touch cuvette anymore

PT:	* 29 mOD	01:48
PT:	* 185 mOD	01:09
PT:		00:00
PT:		01:00

channels are in measurement
Actual optical signal is displayed

PT:	S 31 mOD	01:52
PT:	R 188 mOD	01:13
PT:		00:00
PT:		01:00

a result is found on channel 2 , but it
wait until channel 1 is ready.
Channel 1 is running in high sensitivity

Blinking flag "*" indicates an ongoing measurement

Blinking flag "S" indicates, that the analyzer had switched to the high sensitivity.

Blinking flag "R" indicates, that a result was found , but wait for display and print until all channels are finished.

The instrument will display and results sorted . If the instrument finds a result on optic channel x, it will not display as long as lower numbered channels are in measurement. In this case the actual optical density is flagged with "R" (result found).

5.3.7. Manual break of measurement

To cancel the measurement, press **"Optic"** keys. This will stop the reaction. All optic channels must be not in measurement in order to return to the main menu.

5.3.8. Return to main menu

Press key "MENU" . All optic channels must be not in measurement in order to return to the main menu.
Stop any measurements (refer section above.)

5.3.9. Unit Key Functions

Once the measurement is complete, results are displayed in blinking mode scrolling through all the units (s, %, INR,...) . With a press on any UNIT key the blinking will be stopped on all channels and the units will be scrolled manually.

5.3.10. Stopwatch Functions

Set timer alarm to OFF (See menu SETUP SYSTEM). Press the Timer keys 1-4 to start stopwatch and press again to reset.

5.3.11. Alarm Functions

Set timer alarm to 30s upto 500s (See menu SETUP SYSTEM). Press the Timer keys 1-4 to start countdown and alarm.Press the Timer key again or start channel to reset.

5.3.12. Result Warning Messages

Additionally to the result the instrument may inform the operator at critical samples by attaching status characters or error messages

DISPLAY	PRINTER	MEANING
*	*	Out of calibration (i.e. *167 %)
> <	>	Out of scale (i.e. >999.9 mg/dL)
++++	NO CLOT DETECTED	No clot detected within 300 seconds
----	NO CLOT DETECTED	Clot detected before deadtime
???.?	COAGULATION ERROR	Detected reaction was not valid for coagulation
OPTIC	LOW SIGNAL	Light transmission was not enough.
"*" (blinking)		Ongoing measurement
"S" (blinking)		High sensitivity mode
"R" (blinking)		result found
"B"		Biphasic APTT found (maybe DIC disease)
"F"		Low fibrinogen found (maybe liver disease)


Read section troubleshooting for further details.

5.4. Hidden Functions

To protect the instrument from improper use, critical routines can be called by pressing special key combinations.

5.4.1. Set system to default

Within the "main menu" press key "**OPTIC 1**" and "." and "**ENTER**" simultaneously.



INIT TEST:

- all tests will be correct initiate into the memory. This means in detail that all calibration points will be reset to zero, the method to clotting, the dimensions to seconds.
- The language is set to English
- The format of fibrinogen concentration is mg/dL



This routine is only recommend when a new software will be updated on the instrument

INIT BIOS:

- Integration-time of the receiver to 100 µs
- Temperature Basis Value
- LCD-Contrast to 33
- Speed of Mixer to 215
- Frequency of beep to 325
- The serial interface is configured to 9600 Baud, Stop 1, Data 8, Parity Off



This routine is only recommend when the System pcb (Main board) will be exchanged. The procedure must be followed by a temperature adjustment (see section 5.4.2)!

INIT ALL:

- INIT TEST
- INIT BIOS

5.4.2. Change temperature adjustment


When entering the position "TEMPERATURE" in the submenu "SETUP SYSTEM" the current temperature can be changed by pressing the keys "OPTIC 4" and "UP or DOWN" simultaneously. The temperature will be in- or decrement by 0.1°C. When pressing the keys "OPTIC 4" and "UP or DOWN" and "1" simultaneously, the temperature will be in- or decrement by 1°C. For more information refer to section **8.2 Temperature Adjustment**.

5.4.3. Set all test calibration points to zero

When entering the position "STANDARD CURVE" in the submenu "SETUP TEST" all calibration points can be set to zero by pressing the keys "0" and "ENTER" simultaneously.

5.4.4. Set OD correction

Change to "TEST SETUP" and select test, which you want to correct.
When method is blinking, press key "1" and "MENU" simultaneously.



The measured optical density of the instrument can be corrected by a factor for each test. Therefore it is possible to adapt other reagents.

OD-CORRECTION = 0 no effect (default):

OD-CORRECTION > 0

will reduce sensitivity of method. This is used, if the reagent gives a too high signal response which can cause false early result (eg. PTT is 18s instead of 35s). A too less sensitive setting of the method can cause more "+++" result (no clot found).

Positive OD-CORRECTION < 0 will cause:

will increase sensitivity of method. This is used, if the reagent gives a too low signal response and therefore no clot is detected (+++ results). A much sensitive setting of the method can cause more false early results (eg. PTT is 18s instead of 35s).

OD-CORRECTION = -100%


will invert the signal by -1. This is used for reagents which negative signal response, like a D-dimer assay for red light measured with blue light.



***Improper setting can cause false result. Consult your local distributor or
TECO GmbH before changing the OD correction***

5.4.5. Set COAG correction

Change to "TEST SETUP" and select test, which you want to correct.
When method is blinking , press key "2" and "MENU" simultaneously.



With the Coag Correction the instrument can correct the result for better correlation to other systems or reagents.

New result = factor x result

Factor = 1.00 (COAG CORRECTION = 0%)
Factor = 0.90 (COAG CORRECTION = -10%)
Factor = 1.10 (COAG CORRECTION = +10%)

Example: measured normaltime of PT-S = 15s.

By setting COAG-CORRECTION=-20% the instrument will display 12s.
(New result = 0.80 * 15s = 12s)

Negative COAG-CORRECTION will cause shorter clotting times.

Positive COAG-CORRECTION will cause longer clotting times.



An new calibration run must be followed after every change COAG correction

6.0 Service Menu

1 SYSTEM ANALYSIS
2 ERRORLEVEL
3 OPTIC VALUES

6.1. System Analysis

Press #1 from the Service submenu to enter. Perform a **"SYSTEM ANALYSIS"** to check instrument operational status. In **"SYSTEM ANALYSIS"**, the Coatron M4 check the optic, temperature, Memory values and Analog to Digital Conversion. The Error Level for each channel is determined and displayed if a system error or warning is identified.

An example of a **"SYSTEM ANALYSIS"** screen is:

ERROR IN SYSTEM
OPTIC1: 18 OPTIC2:00
OPTIC3: 00 OPTIC4:00
PRESS ANY KEY

The number after CH1-4 is the error code displayed in hexadecimal format (see 5.2).
Convert the number into binary format and compare bit by bit to the defined errorlevel (see 5.2) .

NOTE: Enter the "SERVICE MENU" and select "ERRORLEVEL" . The errorcodes will then be displayed in binary format. (see 5.2)


Some examples for error codes:

HEX	BINARY	MEANING
00	0000 0000	error free
01	0000 0001	temperature !
08	0000 1000	signal
10	0001 0000	service
40	0100 0000	AD-Conversion !
80	1000 0000	keyboard !
19	0001 1001	temperature + signal + service

After the self-check is complete, a service report is automatically printed. To determine the specific warning or error, enter #2, **"Error Level"**.

Error Level

When the self-check is followed by an error message, a 8 BIT long number in binary code can be seen. Each BIT is indicative of a specific error or warning:



Example: The Error Level "00000001" indicates that the temperature is out of range.

Refer to the troubleshooting guide for information on corrective actions.

6.2. Optic-Values

From the "**SERVICE**" submenu, remove cuvettes and press #3 to enter "**OPTIC VALUES**". The values of signal, noise and service is displayed for each channel.

An example of a "**OPTIC VALUES**" screen is:

CH1:	32432	(202)
CH2:	32283	(191)
CH3:	32169	(213)
CH4:	32574	(207)

Recommended values: *

SIGNAL :	31500 – 33500
SERVICE:	130 – 240




Figure 11 - Flow Diagram of "SERVICE" Submenu

7.0 **Troubleshooting**

TEMPERATURE ERROR	
DESCRIPTION	Temperature is not between 36.0 – 38°C
PROBABLY CAUSE	1. temperature draft environment (window) 2. electronic error 3. temperature is not adjusted correctly
CORRECTIVE ACTION	1. place instrument in a draft free environment without direct sun lighting 2. allow instrument to heat up at least 15 min after temp. adjustment or boot-up 3. adjust temperature (refer sections 8.1 & 8.2) Contact technical services at local distributor if error persists.

OPTIC ERROR : SIGNAL	
DESCRIPTION	This error occurs if not enough light is flashing onto the receiver. The optical signal must be between 20000 - 35000
PROBABLY CAUSE	1. temp. out of range 2. dirt in optics 3. LED's on optic block defective 4. chip errors on boards
CORRECTIVE ACTION	1. wait until instrument is temperatured to 37°C 2. clean optics (refer section 8.3) 3. replace optic block 4. Contact technical services at local distributor if error persists

OPTIC ERROR : SERVICE	
DESCRIPTION	This error occurs if the required signal amplification is not between 110 - 300
PROBABLY CAUSE	1. temp. out of range 2. dirt in optics 3. LED's on optic block defective 4. chip errors on boards
CORRECTIVE ACTION	1. wait until instrument is temperatured to 37°C 2. clean optics (refer section 8.3) 3. replace optic block 4. Contact technical services at local distributor if error persists

SYSTEM ERROR : AD-Conversion	
DESCRIPTION	The photo receiver convert light intensity to analog direct current. This current Is converted to a digital signal. Every digital value has his own signature. If there is lost a value, the system error occurs.
PROBABLY CAUSE	1. CPU defect 2. ADC defect
CORRECTIVE ACTION	Contact technical services at local distributor

ANALYSIS ERROR : + + + +	
DESCRIPTION	No clotting reaction is detected within 300s
PROBABLY CAUSE	1. There is really no clot 2. The fibrinogen concentration in the sample is below 50 mg/dL 3. started channel and pipeted channel are not equal
CORRECTIVE ACTION	1. confirm correct handling 2. Determine clotting time with research software "TECMONI" or observe the optical values on the display during screen. The time until a change in signal can be observed is the clotting time 3. change to recommended reagent provider

ANALYSIS ERROR : ----	
DESCRIPTION	Clotting reaction start and end before deadtime
PROBABLY CAUSE	1. PT based test clot before 7 sec 2. PTT based tests clot before 15 sec
CORRECTIVE ACTION	change to recommended reagent provider

ANALYSIS ERROR : ????	
DESCRIPTION	The instrument detected a reaction, but was not able to verify as a clot reaction.
PROBABLY CAUSE	1. air bubbles 2. touching of cuvette 3. clot reaction starts before deadtime
CORRECTIVE ACTION	1. avoid air bubbles – pipet against cuvette wall 2. avoid touching the cuvette during measurement start 3. change to recommended reagent provider

ANALYSIS ERROR : OPTIC	
DESCRIPTION	The received signal is below 400 digits
PROBABLY CAUSE	1. very turbid samples or reagents 2. dirt in optics 3. optic is defective
CORRECTIVE ACTION	1. clean optics 2. check optic value (refer section 6.2) 3. change to recommended reagent provider

ANALYSIS ERROR : >, <	
DESCRIPTION	Assay range limitation
PROBABLY CAUSE	1. calculated units are not within 0-999 2. calculated activity is not within 3 – 180% 3. calculated ratio is not within 0-9 4. calculated INR is or within 0-12
CORRECTIVE ACTION	Check test setup and recalibrate test if necessary

ANALYSIS ERROR : *	
DESCRIPTION	Calculated result is extrapolated. Extrapolated results can vary to intrapolated results.
PROBABLY CAUSE	The obtained result is outside of calibration
CORRECTIVE ACTION	Extend the calibration curve. High and low standards are helpful.

8.0 **MAINTENANCE**

8.1. Recommended Maintenance

- Daily: Check that the optic unit and the filters are free from dirt. Clean with lint-free tissue paper if necessary.
- Monthly: Check the temperature of the incubator block. When the green LED is on, place water in a reagent bottle and place in one of the reagent positions. Insert a thermometer and record temperature after 10 minutes. The temperature should be in the range of 36.5°C and 37.5°C. Proceed to next section if the temperature is incorrect.
- Yearly: Service check by authorized technical person.
(total cleaning of instrument, function tests, quality-control routines)

8.2. Temperature Adjustment

1. When the green LED is on, fill a reagent container (bottle) with water and place in a reagent position on the incubator block. Place a thermometer in the water.
2. Allow to warm for 10 minutes.
3. Enter the submenu "Setup System" and advance to the temperature screen. The current temperature of the **Coatron M4** is displayed.
4. Compare the temperature displayed by the system and the thermometer. If the temperature is different, adjust the temperature on the **Coatron M4** by simultaneously pressing the Up/Down cursor keys *and* the "OPTIC 4" key. (To increase or decrease using larger increments, press additionally the numeric "1" key.)
5. Wait until a stable temperature of 37.0°C is displayed on the **Coatron M4**. Check and correct the system temperature if not equivalent to the external thermometer.
6. If both the thermometer and instrument display the same temperature, press the "Enter" key and exit back to the main menu.

8.3. Cleaning procedures

The casing should be cleaned with a light alcoholic solution or with a light soap solution, using a soft foam only. The optical filter should be cleaned with just light alcoholic solution. Dry well with clean, fat-free wipes

Warning:

- ***Do not use aggressive cleaning solutions !***
- ***Do not use hard foam material or metallic foams.***

9.0 APPLICATIONS

A summary of the current test applications supported by the **Coatron M4**, the choice of medical dimensions and the numeric test codes are shown below:

Test	medical dimensions	
	Unit	# CODE
STAT	s - % - R - I - U	00
PT +/- FIB	s - % - R - I - E - U	01
PTT	s - R	02
TT	s - R	03
FIB(CLAUSS)	s - U	04
VT	s - R	05
F2	s - %	06
F5	s - %	07
F7	s - %	08
F8	s - %	09
F9	s - %	10
F10	s - %	11
F11	s - %	12
F12	s - %	13
F15 (Fletcher)	s - %	14
HEP	s - U	15
AT3	s - % - U	16
APCR	s - R	17
PS	s - % - U	18
PC	s - % - U	19
LA-S	s - R	20
LA-C	s - R	21
DD	E - U	22
VWF	E - U	23
ECAH	E - % - U	24
ECAT	E - % - U	25
PLG	E - % - U	26
a2AP	E - % - U	27

GLOSSARY

- ◆ E, result in optical density (Extinction); precision XXXX
- ◆ s, result in seconds; precision XXX.X
- ◆ %, result in activity; precision XXX.X
- ◆ U, result in mg/dl, mg/l, µg/ml, g/l; precision XXX (except HEP: X.XXX)
- ◆ R, result in ratio; precision XX.XX
- ◆ I, result in INR; precision XX.XX

All application can be run with a minimum volume of 75 µL. But in some cases it is recommend to use higher volume levels to keep the precision in a high accuracy. The next pages suggest test procedures, which are tested on the instrument and which guaranty a Coefficient Variation (CV) of below 5 % (depend on test, pipette).

9.1. Prothrombin Time

REAGENT PREPARATION

Reconstitute Thromboplastin reagents according to the package insert.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check setup system if necessary.
5. Check setup test if necessary . Enter new calibration curve data to obtain results in % activity; enter ISI of thromboplastin reagent and in-house determined mean normal PT if INR results are desired. Enter in-house determined mean normal value if results in R are required.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #01. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

Clotting Method:

1. Pipette **25 μ L plasma** to cuvette(s).
2. Prewarm plasma for 2 min, or the time indicated in the reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
3. Place Thromboplastin reagent with stir bar in large central reagent position.
4. Transfer cuvette to measuring position.
5. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.
6. Add **50 μ L prewarmed Thromboplastin reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test). Repeat for remaining channels.

7. The instrument will read for 300 secs. If no clot is detected, the display will read "***".
8. The result is displayed in seconds. Press the corresponding "Unit" key for conversion of results.

ASSAY CALIBRATION

For calibration curves, a minimum of two values is required, with a maximum of 5. **It is highly recommended that more than two calibration points be used.**

1. Make dilutions of Normal Coagulation Reference Plasma I in 0.9% NaCl (saline). Refer to the table below for preparation of standards
2. Determine PT with undiluted normal plasma in duplicate.
3. Determine PT with diluted plasmas in duplicate.
4. If no clot time is detected for the diluted plasma samples (<25%), change the sensitivity in "Setup Test" to "High" and repeat.
5. Enter calibration data, % activity and seconds.
6. Check calibration curve with different controls, Controls Level I, II & III.

% Activity	Dilution	Preparation
100	none	
50	1:2	200 uL plasma + 200 uL saline
25	1:4	100 uL plasma + 300 uL saline
12.5	1:8	50 uL plasma + 350 uL saline

9.2. Derived Fibrinogen

REAGENT PREPARATION

Reconstitute Thromboplastin reagents according to the package insert.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary
5. Check test setup if necessary. Select method clotting+fibrinogen, enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #1. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

Clotting method with fibrinogen:

1. Pipette **25 μ L plasma** to cuvette.
2. Prewarm plasma for 2 min, or the time indicated in reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
3. Place Thromboplastin reagent with stir bar in large central reagent position.
4. Transfer cuvette to measuring position.
5. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.
6. Add **50 μ L prewarmed Thromboplastin reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test). Repeat for remaining channels.
7. The instrument will read for 300 secs. If no clot is detected, the display will read "****".
8. The result is displayed in seconds. Press the corresponding "Unit" key for conversion of results.

ASSAY CALIBRATION

For the derived fibrinogen test a calibration curve must be entered. **It is highly recommended that more than two calibration points be used.**

1. Test fibrinogen reference plasmas in duplicate. For example, test fibrinogen reference plasmas with low (~100 mg/dL), normal (~250 mg/dL) and high (>350 mg/dL) fibrinogen concentration. Record the Extinction (OD) values for all samples.
2. Enter the fibrinogen data for the calibration curve, Extinction & fibrinogen values.
3. Check the calibration curve with controls.
4. Note - If % PT *and* derived fibrinogen values are required, two calibration curves must be entered.

9.3. Clauss Fibrinogen Assay

REAGENT PREPARATION

Reconstitute Thrombin reagent 100 NIH U/mL according to the package insert.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary.
5. Check test setup if necessary and enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select Fib with the Up/Down arrow keys or enter the numeric test code, #6. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

All quality control and patient samples are diluted 1:10 in Imidazole Buffered Saline (IBS) for testing. If the clotting times fall outside of the linear curve, prepare and test 1:5 or 1:20 dilutions as needed.

1. Pipette **50 μ L of diluted sample** to cuvette.
2. Prewarm sample for 5 min, or the time indicated in reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
3. Transfer cuvette to measuring position.
4. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.
5. Add **25 μ L 100 NIH U/mL Thrombin reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test).
6. The instrument will read for 300 secs. If no clot is detected, the display will read "+++" and "No Clot Detected" will print.
7. The result is displayed in seconds, and both this result and the fibrinogen concentration are automatically printed. press the corresponding "Unit" key for fibrinogen concentration if a printer is not attached.

8. For samples diluted 1:10, this is the final result. For other dilutions, the result must be corrected. For example, if the sample was diluted 1:5, divide the result by 2; if the sample was diluted 1:20 or 1:40, multiply the result by 2 or 4, respectively.

ASSAY CALIBRATION

1. Prepare standards using Normal Coagulation Reference Plasma. A *suggested* standard curve is shown in the following chart:

Dilution	Preparation
1:5	200 μ L plasma + 800 μ L IBS
1:10	100 μ L plasma + 900 μ L IBS
1:20	50 μ L plasma + 950 μ L IBS
1:35	100 μ L plasma + 3.4 mL IBS

2. Mix the standards and assay each in quadruplicate.
3. Enter the calibration data in the "Setup Test" submenu. For example: the assigned fibrinogen value for the plasma used to prepare the standard curve is 250 mg/dL. The 1:10 dilution corresponds to 100% activity, therefore it is equal to 250 mg/dL. The 1:5 dilution is twice as concentrated, it is equal to 500 mg/dL. (The 1:20 dilution = 125 mg/dL and the 1:35 = 71.4 mg/dL.)
4. Verify the calibration curve with different Control Plasmas Level I,II,III.

9.4. Thrombin Time Assay

REAGENT PREPARATION

Reconstitute Thrombin Reagent 10NIH/vial according to the package insert. (use "europe" procedure with concentration of 5NIH/ml)

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary
5. Check test setup if necessary. Select method clotting+fibrinogen, enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #1. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

Clotting Method:

1. Pipette **50 μ L of undiluted plasma** to cuvette.
2. Prewarm sample for 3 min, or the time indicated in reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
3. Transfer cuvette to measuring position.
4. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.
5. Add **50 μ L Thrombin Time Reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test).
6. The instrument will read for 300 secs. If no clot is detected, the display will read "+++" and "No Clot Detected" will print.
7. The result is displayed in seconds and is automatically printed.

9.5. APTT

REAGENT PREPARATION

Refer to package insert for APTT and CaCl_2 reagents.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary
5. Check test setup if necessary. Select method clotting+fibrinogen, enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #1. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

Clotting Method:

1. Pipette **25 μL plasma** to cuvette.
2. Prewarm plasma for 2 min, or the time indicated in reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
3. Place CaCl_2 in large central reagent position.
4. Add **25 μL of the APTT reagent**. Incubate for 3 or 5 minutes (refer to package insert).
5. Transfer cuvette to measuring position.
6. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.
7. Add **25 μL prewarmed CaCl_2 reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test).
8. The instrument will read for 300 secs. If no clot is detected, the display will read "****".
9. The result is displayed in seconds. To obtain results in R, press the "Unit 1" key.

9.6. PT-Based Factor Assays (II, V, VII & X)

REAGENT PREPARATION

Reconstitute Factor Deficient Plasma and Thromboplastin reagent according to package insert.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary
5. Check test setup if necessary. Select method clotting+fibrinogen, enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #1. If any warning or error message appears, refer to section 7.0.

TEST PROCEDURE

All quality control and patient samples are diluted 1:10 in Imidazole Buffered Saline (IBS) for testing. If the clotting times fall outside of the linear curve, prepare and test 1:5 or 1:20 dilutions as needed. It is recommended to test 2 dilutions (1:10 & 1:20) for each sample. Refer to the reagent package insert for more information.

Clotting Method

1. Place Thromboplastin reagent with stir bar in large central reagent position.
2. Pipette **25 µL of diluted plasma and 25 µL of deficient plasma** to each cuvette. Refer to chart for sample preparation.
3. Incubate for 2 min, or the time indicated in reagent package insert. Press the "TIMER 1" key to start stop-watch 1.
4. Transfer cuvette to measuring position.
5. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing "OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.

6. Add **50 µL prewarmed Thromboplastin reagent** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test). Repeat for remaining channels.
7. The instrument will read for 300 secs. If no clot is detected, the display will read "+++" and "No Clot Detected" will print.
8. The result is displayed in seconds. Press the corresponding "Unit" key for conversion of results if a printer is not attached.
9. For patient and control samples diluted 1:10, this is the final result. If other dilutions are tested, the calculated value should be multiplied by the appropriate dilution correction factor. (i.e., samples diluted 1:20, multiply result by 2; for 1:40 dilutions, multiply by 4, etc.)

ASSAY CALIBRATION

For calibration curves, a minimum of two values is required, with a maximum of 5. **It is highly recommended that more than two calibration points be used.**

1. Make dilutions of Normal Coagulation Reference Plasma in Imidazole Buffered Saline (IBS). A *suggested* standard curve is shown below.
2. Assay standards in quadruplicate as described.
3. Enter calibration data (% activity and seconds) in "Setup Test". Check calibration curve with different controls.

Sample	Dilution	Preparation
100% Standard	1:10	100 uL reference plasma + 900 uL IBS
50% Standard	1:20	50 uL of reference plasma + 950 uL IBS
25% Standard	1:40	25 uL of reference plasma + 975 uL IBS
12.5% Standard	1:80	500 uL of 25% standard + 500 uL IBS
Patient or Control	1:10	100 uL sample + 900 uL IBS

9.7. APTT-Based Factor Assays (VIII, IX, XI & XII)

REAGENT PREPARATION

Reconstitute Factor Deficient Plasma and APTT reagents according to package insert.

SYSTEM PREPARATION

1. Turn on instrument and wait for green LED light to come on.
2. Turn on printer if connected.
3. Connect optional Autopipette to system.
4. Check system setup if necessary
5. Check test setup if necessary. Select method clotting+fibrinogen, enter new calibration curve data.
6. Return to main menu and enter "Analysis" by pressing #1. Select PT with the Up/Down arrow keys or enter the numeric test code, #1. If any warning or error message appears, refer to section 7.0..

TEST PROCEDURE

All quality control and patient samples are diluted 1:10 in Imidazole Buffered Saline (IBS) for testing. If the clotting times fall outside of the linear curve, prepare and test 1:5 or 1:20 dilutions as needed. It is recommended to test 2 dilutions (1:10 & 1:20) for each sample. Refer to the reagent package insert for more information.

Clotting Method

1. Place CaCl_2 in large central reagent position.
2. Pipette **25 μL of diluted plasma and 25 μL of deficient plasma** to each cuvette. Refer to chart for sample preparation.
3. Incubate for 2 min, or the time indicated in reagent package insert.
4. Add **25 μL of APTT reagent**. Refer to package insert for appropriate activation times. Press the "TIMER 1" key to start stopwatch 1.
5. Transfer cuvette to measuring position.
6. While incubating, press "OPTIC 1". If selected, enter PAT-ID with numeric keys or Up/Down keys. Confirm by pressing

"OPTIC 1" again. The message "ACTIVE" is displayed and channel 1 is ready to start the reaction. Repeat for the remaining channels.

7. Add **25 μL prewarmed CaCl_2** and simultaneously press the "OPTIC 1" key. The test will automatically start if using the Autopipette. (CAUTION: When the test procedure is running, pressing the "OPTIC 1" and the "Enter" keys will interrupt the test). Repeat for remaining channels.
8. The instrument will read for 300 secs. If no clot is detected, the display will read "+++" and "No Clot Detected" will print.
9. The result is displayed in seconds. Press the corresponding "Unit" key for conversion of results if a printer is not attached.
10. For patient and control samples diluted 1:10, this is the final result. If other dilutions are tested, the calculated value should be multiplied by the appropriate dilution correction factor. (i.e., samples diluted 1:20, multiply result by 2; for 1:40 dilutions, multiply by 4, etc.)

ASSAY CALIBRATION

For calibration curves, a minimum of two values is required, with a maximum of 5. **It is highly recommended that more than two calibration points be used.**

1. Make dilutions of Normal Coagulation Reference Plasma in Imidazole Buffered Saline (IBS). A *suggested* standard curve is shown below.
2. Assay standards in quadruplicate as described.
3. Enter calibration data (% activity and seconds) in "Setup Test". Check calibration curve with different controls.

Sample	Dilution	Preparation
100% Standard	1:10	100 uL reference plasma + 900 uL IBS
50% Standard	1:20	50 uL of reference plasma + 950 uL IBS
25% Standard	1:40	25 uL of reference plasma + 975 uL IBS
12.5% Standard	1:80	500 uL of 25% standard + 500 uL IBS
Patient or Control	1:10	100 uL sample + 900 uL IBS

10.0 **SPECIAL FUNCTIONS**

10.1. Software Upgrading


The process of evolution never stops. This is also valid for technical systems. Complaints or ideas from the customers and their direct solution or realization force an instrument to more and more functionality and stability. The **Coatron M4** supports a very easy procedure to use a new revision of software. For more information about latest software revision, ask your local distributor .

Requirements :

- Windows PC ,
 - download cable ,
 - update files - flash.bat , download.exe , flashdl.exe , flashv25.exe and the analyser firmware xxxxx.img (i.e. c4u0111.img → Coatron M4 1.11 400nm)
1. Save any important test setups of the instrument to paper
 2. Reboot your PC. This will reset your comm ports !
 3. Extract „Flash Disk to any local folder on your PC.
 4. Connect Coatron to PC with 0 Modem cable (PIN 2-3-5 to PIN 3-2-5)
 5. Double click to "upgrade.bat" and follow dialogue.
 6. After a few minutes the Coatron will reboot.
 7. Re-Initialisize the memory
 - a. Power off/on instrument
 - b. Press key "**Optic 1**" + "." + "**ENTER**" **simultaneously** -> message "INIT SYSTEM..." should appear and afterwards instrument should reboot automatically
 8. Restore test setup parameter manually

Problems during updates:

-Update stays at following screen



```
C:\WINNT\system32\cmd.exe
D:\TECO\BAUMI\CA4\PROGRAMM>Download -c1 -iflashCA4 -y -q
TECO-DOWNLOAD TestVersion 2.00 Sep 14 1998
(C) DOLI Elektronik GmbH.(1995-1998)

MinimalDownload
Start-Segment (Hex) : 100
Serieller Download <115200 Baud>
Serieller Download <DE-Adr=00000009>

D:\TECO\BAUMI\CA4\PROGRAMM>d1CA4 -c1 -iCA4.bin -eC000 -q $r $x $v

DOLI FLASH-Programmierung
FLASHDL Version 2.03 Jun 2 1999
(C) DOLI Elektronik GmbH.(1995-1999)

Infile CA4.bin
Backupfile flash.sav
Kommunikation COM1, 115200 Baud
Nur EraseBlock sichern!!!
NoSave!!!
NoVerify!!!

!tarte Sync. zum SUB.
```

- ➔ Reboot PC and try again
- ➔ If transfer still failed, boot PC with MS DOS floppy disk and perform MS DOS update

- Nothing happened after click to upgrade.bat
 - ➔ perform update with MS DOS commands

- MS DOS update (for Windows 2000 / XP)
 - open data explorer and extract all files maybe to drive "C" , folder TECO
 - goto MS Windows "START/RUN" and enter "cmd" -> DOS command window will open
 - change to directory with command "cd c:\teco"
 - enter command "dir" -> all files including upgrade.bat should be listed
 - enter command "upgrade.bat" -> a dialog should open.
 - Follow dialog. Update will start and when finished the instrument will reboot.

11.0 **TECAM SMART the LIMS solution**

11.1. General

TECAM SMART is a user friendly software with LIMS functionality (Laboratory Information Management. It allows to sample all instrument results and manage it in a database application. The software can be downloaded from TECO homepage www.teco-gmbh.com . The trial version allows only saving 5 results each start but has no other limits. With a purchased license it can be switch to a registered version.

Requirements:

- Pentium III , 128MB Ram
- Microsoft Windows 2000 or XP

DataBase:

- Microsoft Access Jet Engine 4.0
- Max. 2GB (200.000 results)

Supported instruments:

- Teco Coatron M1 V1.18 or higher
- Teco Coatron M2 V1.13 or higher
- Teco Coatron M4 V1.08 or higher
- Teco Techrom IV V1.20 or higher

Key features:

- Account Management: operator can login as Administrator or user in order to control and manage the database with different access rights
- DataBase Management: operator can create, backup, import or export database
- Reaction Curve: can be displayed in fix or autoscale mode
- Patient DataBase: PID can be linked to patient names
- Filter functions: Every record field can be filtered (e.g. last 30 days all PT)
- Report Generator: Day report , patient report , or last x days report can created just by filter the fields
- Statistical analysis: includes mean, SD,CV for filtered data
- Online Update: software be updated by internet

11.2. Interface Protocol

Protocol: TECAM = TECO Coatron Monitoring
Instruments: Coatron M2 V1.13b or higher
Coatron M4 V1.08b or higher
Techrom IV V1.20 or higher
Interface: 9600 Baud , no parity , 8 bit . 1 stop bit
Transmission: result data = yes
Reaction curve data = yes

ASCII command characters:

STX	start of transmission	asc(2)
ETX	end of transmission	asc(3)
TAB	vertical tabulator	asc(9)
LF	line feed	asc(10)
CR	carriage return	asc(13)
DLE	idle	asc(16)

Overview of records and principle:




Figure 12 – LIMS communication

The result record is always transmitted from the Coatron.
The curve record transmission must be activated from PC with a request record.

The curve request record:

Function: enables curve transmission !
Definition: <DLE &"4" CR>

The curve undo record:

Function: disables curve transmission !
Definition: <DLE &"3" CR>

The curve start record

Function: Indicates that a new measurement had been started.
The transmission need to be enabled by the "Curve Request Record".

Definition: <STX &"P"><Optic><LF & ETX>

Fields: Delimiter: asc(9) = tabulator
<STX & "P"> Record header ; STX=asc(2)
<Optic> optic cannel 1-4 (channel 1 = 0)
<LF & ETX> Record end; LF=asc(10) , ETX=asc(3)

The curve record

Function: Send actual optical density every second
The transmission need to be enabled by the "Curve Request Record".

Definition: <STX &"O"><Optic><mOD ><LF & ETX>

Fields: Delimiter: asc(9) = tabulator
<STX & "O"> Record header ; STX=asc(2)
<Optic> optic cannel 1-4 (channel 1 = 0)
<mOD> milli optical density
<LF & ETX> Record end; LF=asc(10) , ETX=asc(3)

Example: Signal = 0,999 E in channel 1--> <STX "O"><0><999><LF & ETX>

The result record

Function: send after every result

Definition: <STX & "R"> <Optic> <SID> <Error> <Date> <Flags> <PID> <Test>
<sec> <mOD> <%> <R> <INR> <Unit> <Scale> <LF & ETX>

Fields:	Delimiter:	asc(9) = tabulator
	<STX & "R">	Record header ; STX=asc(2)
	<Optic>	optic cannal 1-4 (channel 1 = 0)
	<SID>	system identifier
	<Error>	ystem error , hexadecimal
	<Date>	ddmmyy or 000000 if not used
	<Flags>	Result flags : + , - , * , >
	<PID>	Result or patient ID number
	<Test>	name of test ("PT", "APTT", ...)
	<sec>	125 (no decimal point. Divide by 10 -> 12.5s)
	<mOD>	000 (no deciaml point)
	<%>	100.0 (with decimal point)
	<R>	1.00 (with decimal point)
	<INR>	1.00 (with decimal point)
	<Unit>	300. (with decimal point)
	<Scale>	mg/dL , ng/mL ,
	<LF & ETX>	Record end; LF=asc(10) , ETX=asc(3)

Examples:

Test PT with 12.5s 1.02INR 98.5% PID=1000 ERR=0 Kanal=1
<STX R><0><1011000><0><311204><0><1000><PT>
<125><0><98.5><1.02><0><><LF ETX>

Test FIB with 10.0s 200 mg/dL PID=1000 ERR=0 Kanal=1
<STX R><0><1011000><0><311204><0><1000><FIB>
<100><0><0><0><200><mg/dL><LF TX>

Test AT3with 0.651E 78% PID=1000 ERR=0 Kanal=1 (chromogen)
<STX R><0><1011000><0><311204><0><1000><AT3>
<0><651><78><0><0><><LF ETX>

12.0 **PRODUCT CATALOGUE**

Coatron M4 (with Standard Package)

Cat. Number: 80 420 00

including:

- 1 pc Power supply (93-240Vac/50-60Hz)
- 1 pc Mains cable
- 25 pcs Double cuvettes (2 wells/each)
- 5 pcs Reagent container
- 1 pc Stirring magnets
- 1 pc Reagent adaptor (Ø22.5mm)
- 1 pc Printer cable
- 1 pc Pipette fastener (with 2 screws)
- 1 pc Operation manual
- 1 pc Warranty card
- 1 pc Installation report
- 1 pc Service report

Consumables and Accessories

Cat.No.	Description	Content	Qty
19 000 02	Double-cuvette (2 pos/ea)	250	1 Pac
19 029 00	4-stage-Autopipette 25/50/100/200 µl	1	1 Pc
19 030 01	Reagent adaptor Øi=22,5mm	1	1 Pc
19 030 02	Reagent adaptor Øi=22,8mm	1	1 Pc
19 030 03	Reagent adaptor Øi=24,2mm	1	1 Pc
19 030 04	Reagent adaptor Øi=27,8mm	1	1 Pc
19 030 05	Reagent adaptor Øi=25,2mm	1	1 Pc
19 030 06	Reagent adaptor Øi=18,6mm	1	1 Pc
19 070 00	CCD Barcode Reader	1	1 Pc
80 050 06	Stirring magnets	4	1 Pac
80 050 09E	Power cord, Europe	1	1 Pc
80 050 09U	Power cord, US	1	1 Pc
80 525 00	Reagent container Ø 22,5mm	100	1 Pac
80 530 00	Reagent tubes Ø 16mm	100	1 Pac
80 535 00	Thermal printer (EU)	1	1 Pc
80 920 07	Pipette tips 25/50/100/200 µl	1.000	1 Pac
80 920 12	Thermal paper for printer	5	1 Pac
80 920 14	Printer cable for printer	1	1 Pc
80 997 00	TECAM Smart Software Patient Management, Monitoring Research CD and RS232 cable	1	1 Set