IVD

REF

A0260-050

Intended Use

Clotting test for the quantitative determination of prothrombin time (PT) in citrated human plasma or whole blood according to Owren. Intended to be used by professional laboratory personnel using coagulation analysers. Is not intended for self testing/POCT. The determination of PT is used for the global evaluation of the extrinsic pathway of coagulation.

Prolonged clotting times may be observed in the following situations: deficiency or inhibition of coagulation factors II, VII, X, Vitamin-K, presence of Vitamin-K inhibitors like Phenprocoumon, Warfarin or others.

Contents & Determinations

Product	TEClot PT-B
REF	A0260-050
TEClot PT-B reagent	5 x 10mL
CaCl ₂ (25mM)	5 x 5mL
Diluent	5 x 5mL

Determinations

Coatron ECO	333 Det.
Amelung KC	250 Det.
Owren's	500 Det.

Ingredients

TEClot PT-B reagent	Thromboplastin from rabbit brain and bovine plasma fraction. The plasma is depleted from coagulation factors II, VII and X, but is rich in fibrinogen and coagulation factor V.
CaCl ₂	25mM calcium chloride solution with 4mM sodium acid (NaN ₃) as preservative
Diluent	High purified water with 4mM NaN ₃

Preparation

TEClot PT-B must be reconstituted in two phases to avoid flakes in the vial.

Phase-1	- Reconstitute TEClot PT-B with Diluent - Let the reagent dissolve for 10 minutes before adding CaCl ₂
Phase-2	- 10min after phase-1 add Calcium chloride

Wait at least 60min until first measurement. Keep reagent at room temperature and mix gently every 15min.

Storage & Stability

Unopened reagents are stable until the expiration date shown on the label stored at 2°-8°C. Reconstituted reagent & control must be stored always in closed condition to avoid evaporation effects. Avoid direct sunlight during storage:

After reconstitution	2-8 °C	20-25 °C	37°C
PT-B Reagent	5 days	5 days	8 hours
Diluent/CaCl ₂	Ready to use. Expiry date see on label.		

Specimen collection and storage ¹

Citrated whole blood:

- Obtain venous blood by clean vein puncture.
- Immediately mix 9 parts blood with 1 part 3.2% sodium citrate (0.105M).
- Analyse within 4 hours.

Citrated plasma:

Centrifuge Citrated whole blood at 1500g for 10 min. (platelet < 10000/μL)

Precautions

- Avoid contact with skin and eyes. Wear suitable protective clothing. All components are checked for HIV, HBV and HCV. However, products from human blood should be considered as potentially infectious
- Do not drain CaCl₂ and diluent into systems with lead or copper. The containing NaN₃ will react to these metals to form explosive compounds

Disposal

Dispose components in accordance with local regulation and refer to the product safety declaration for the link to appropriate hazard and precautionary statements.

Warranty

This product is warranted to perform in accordance with its labelling and literature. TECO disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event will TECO be liable for any consequential damages arising out of aforesaid express warranty

Coatron X ECO procedure (not for other Coatron versions !)

- Pipette 150 μL "PT-B Reagent" into empty cuvette. Incubate 1-5 min at 37°C.

Symbols key:

	Expiry date		In Vitro Diagnostica		Biological hazard		Catalogue Number		Reconstitute with dest. water		Consult accompanying documents
	Store at 2-8°C		EU conformity		Manufacturer		Lot. Number		Ready to use		Authorized Representative

- Add 15 μL "whole blood" (citrate or capillary) or control and immediately mix 10-15x with pipette. Avoid bubbles and remove pipette in pressed condition.

Preparation of calibration:

- 100% Standard: Reconstitute TECal N with 1.7mL
- 25% Standard: Dilute 100μL TECal N with 500μL Owren buffer

Amelung KC procedure

- Add 100μL Owren's buffer and steel ball into empty cuvette and incubate 60sec
- Pipet 20 μL whole blood (citrate or capillary) or TEControl
- Add 200μL PT-B Reagent, preheated to 37°C

Preparation of calibration:

- 100% Standard: Reconstitute TECal N with 1.4mL
- 50% Standard: Dilute 100μL TECal N with 200μL Owren buffer
- 25% Standard: Dilute 100μL TECal N with 500μL Owren buffer
- 12.5% Standard: Dilute 100μL TECal N with 1100μL Owren buffer

Owren's procedure (plasma only !)

- Predilute plasma 1:7 in Owren's buffer (100μL plasma + 600μL buffer)
- Pipet 50 μL prediluted plasma into empty cuvette and incubate 60s
- 100 μL PT-B Reagent, preheated to 37°C

Preparation of calibration:

- 100% Standard: Reconstitute TECal N with 1.0mL
- 50% Standard: Dilute 100μL TECal N with 100μL Owren buffer
- 25% Standard: Dilute 100μL TECal N with 300μL Owren buffer

Expected Results

Normal values: 18 – 26 sec 70 - 130% 0.85 – 1.15 INR
Therapeutic range: 2.00 – 3.00 INR

However, results may vary by instruments, technique, calibration, etc. Each laboratory is recommended to establish its own range on the specific instrument used.

Quality Control

Dimex procedure: Reconstitute TEControl with 1.7mL dest. water
KC procedure: Reconstitute TEControl with 1.4mL dest. water
Owren's procedure: Reconstitute TEControl with 1.0mL dest. water

TEControl can be frozen one time after reconstitution. 120-150 μL stored in closed polypropylen tubes at -20°C is stable for 30 days. TEControl or other commercial control plasma should be used for reliable quality control of performance at a frequency in accordance with good laboratory practice (GLP).

Calibration

INR-System: Run normal control and transfer clotting time as normal time into instrument. Enter correct ISI value into instrument. ISI value is stated on the PT-B vial.

%-System: Run calibrators (100%, 50%, 25%, ...) and transfer results into calibration curve of instrument. (eg. 100% = 22s; 50%=35s; 25%=60s; etc.)
(Verify accuracy of calibration with control plasma normal and abnormal.)

Performance Characteristics**Typical performance on instrument Dimex (whole blood)**

Precision: CV% (within run) CV% (inter-runs)
Normal control < 7,5 < 10,0
Abnormal control < 7,5 < 10,0
(Precision depends very much on pipetting techniques - especially at mixing)

Limitations and interfering substances

- Use only plastic tubes or siliconised glass
- Reagent must be free of flocks. Check for precipitation before use and in case follow note-2.
- Dissolved reagent should be always stored in closed condition at 2-8° and gently mixed before use. Best is to transfer only a certain amount of reagent in to a separate closable cup, which can be kept at room temperature for one day.
- Do not re-use prewarmed reagent on the next day.
- Use only fresh tips to avoid any contamination for the reagent.
- Incubate cuvettes max 5min to avoid evaporation effects
- Check pipette regularly for accuracy
- Check reagent visually for precipitation before use and follow point-2 in see chapter
- Use only samples with less than 1.0 U/mL heparin, 40mg/dL bilirubin, 1000 mg/dL hemoglobin, 1000 mg/dL triglycerides

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