IVD

REF

A0320-050

Intended Use

Clotting test for quantitative determination of the Activated Partial Thromboplastin Time (APTT) in citrated human plasma using silicate as contact activator for factor XII. Intended to be used by professional laboratory personnel using coagulation analysers. The determination of the APTT is used for the global evaluation of the intrinsic pathway and detecting deficiencies of the intrinsic coagulation factors VIII, IX, XI, XII, and Fletcher Factor or other coagulation methods where an APTT reagent is required^{1,2}.

The APTT reagent in the kit contains phospholipids and silica to ensure a highly consistent and stable product³. The APTT reagent is lupus anticoagulant insensitive. Lupus anticoagulant insensitive reagents yield more reliable factor assay results than reagents, which are sensitive to lupus inhibitors⁴.

Prolonged clotting times may be observed in the following situations: deficiency of intrinsic coagulation factors, presence of heparin or other anticoagulants, which affect the intrinsic pathway and in liver diseases.

Contents

Product	TECLOT aPTT-S
REF	A0320-050
aPTT-S reagent	10x5 mL
Determinations*	2000

*Micro method (75µL in total)

APTT-S reagent contains colloidal silicate with phospholipids, buffer and preservatives.

Recommended additional material (not included in package)

Auxiliary reagents	A0350-050 Calcium Chloride 0.025M, 10 x 5mL
Calibration	not required
Quality Control	P6001-010 TECO Control N, 10 x 1mL P6101-010 TECO Control A, 10 x 1mL

Preparation

Ready to use. Swirl APTT reagent gently prior usage

Storage & Stability

Unopened reagents are stable until the expiration date shown on the label stored at 2°-8°C. Opened reagent:

	2-8 °C	20-25 °C	37 °C
APTT-S Reagent	30 days	8 days	8 hours

Precautions

The reagent contains sodium azide (less than 0.1%) to prevent microbial growth. Avoid contact with skin and eyes. Wear suitable protective clothing. Dispose components in compliance with local regulations for infectious material.

Specimen collection and storage⁵

1. Obtain venous blood by clean vein puncture.
2. Immediately mix 9 parts blood with 1 part 3.2% sodium citrate (0.105M) and mix well
3. Centrifuge the specimen at 1500g for 10 min. (platelet < 10000/µL)
4. Separate plasma after centrifugation and store in plastic or siliconised glass tube.
5. Use plasma within 4 hours, otherwise store frozen and thaw just prior to use.

Stability of plasma: 4h at 18-26°C 8h at 2-8° 14d at -20°C 6m at -70°C

Procedure**A. Automated Method: Coatron A**

See application book of device

B. Manual Method: Coatron X

1. Prewarm **CaCl₂** (0.025M) at 37°C for at least 10 min
2. Pipette **25 µL of sample** into a test cuvette. Prewarm at 37°C for 1-2 minutes.
3. Add **25 µL APTT-S reagent** and incubate exactly for **3 min** at 37°C.
4. Add **25 µL of CaCl₂** (0.025M) and simultaneously start test.
5. Record the clotting time in seconds.

Expected Results

Typical normal results are 27-42 sec. However results are influenced by the method of clot detection and can vary from laboratory to laboratory. Each laboratory is recommended to establish its own normal range on the specific instrument used.

Quality Control

TECO Control or other commercial control plasma should be used for reliable quality control of performance at a frequency in accordance with good laboratory practice (GLP). TECO Control can be frozen one time after reconstitution. 120-150 µL stored in closed polypropylen tubes at -20°C is stable for 30 days

Limitations**A. Specimen Collection. AVOID:**

1. Use only plastic tubes or siliconised glass.
2. Delayed mixing of blood with anticoagulant.
3. Contamination with tissue thromboplastin.
4. Improper ratio of anticoagulant with blood.
5. Hemolyzed, icteric or lipemic samples may interfere optical systems

B. Laboratory Techniques

1. Perform tests at 37°C.
2. Use only high purity water.
3. Optimum pH is 7.0-7.5.

C. Interfering substances

1. Bilirubin 40mg/dL
2. Haemoglobin 1000 mg/dL

Performance Characteristics**Typical performance on instrument Coatron X**

Precision: CV% (within run) CV% (inter-runs)

QC control < 3,0 < 5,0

Factor & Heparin sensitivity:

Factor (%)	APTT Clotting time (s)		
	F VIII	F IX	F XI
< 1%	100	80	103
10%	53	52	58
40%	40	39	41
100%	35	35	35

Heparin (U/mL)	0 U/mL	0,2 U/mL	0,4 U/mL
APTT clotting time (s)	35	70	180

These values should be used as guidelines only. Each laboratory should establish factor or heparin sensitivity using its own instruments and techniques.

Warranty

This product is warranted to perform in accordance with its labelling and literature. TECO disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event will TECO be liable for any consequential damages arising out of aforesaid express warranty.

References

1. Proctor RR, Rapaport SI. The partial thromboplastin time with kaolin. A simple screening test for first stage plasma clotting factor deficiencies. *Am J Clin Pathol* 36, 212-219 (1961).
2. Triplett DA, Harns CS, Koepke JA. The effect of heparin on the activated partial thromboplastin time. *Am J Clin Pathol* 70, 556-569 (1978).
3. Stevenson KJ, Easton AC, Thomson JM, Poller L. The reliability of activated partial thromboplastin time methods and the relationship to lipid composition and ultrastructure. *Thromb Haemost* 55, 250-258 (1986).
4. Denis-Magdelaine A, Flahault A, Verdy E. Sensitivity of sixteen APTT reagents for the presence of lupus anticoagulants. *Haemostasis* 25, 98-105 (1995).
5. NCCLS: Guidelines for the Standardized Collection, Transport and Preparation of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays

Symbols key:

	Expiry date		In Vitro Diagnostica		Biological hazard		Catalogue Number		Consult accompanying documents
	Store at 2-8°C		EU conformity		Manufacturer		Lot. Number		Authorized Representative

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