Factor II Deficient Plasma

TECO

In-Vitro-Diagnostic reagent for professional use

IVD

REF

P5001-010

Intended Use

This product is intended for the quantitative determination of factor II in human plasma by clotting assay. Factor II or prothrombin, is a vitamin K-dependent clotting factor which is the inactive precursor of thrombin. Thrombin is generated from prothrombin upon activation by factor Xa. This reaction takes place on a phospholipid surface in the presence of calcium ions and requires factor V as a cofactor. Samples with factor II deficiency have prolonged prothrombin time (PT).

Contents & Determinations

Cat.No.	P5001-010
Factor II Deficient Plasma*	10x1 mL

Determinations

Coatrom M*	Coatrom A4	Coatrom A6
400	200	250

*contains freeze-dried factor II depleted human plasma

**Micro method (75µL in total)

Preparation

Reconstitute each vial with 1 mL distilled water. Mix gently and allow standing at room temperature for 20 minutes for complete dissolution.

Further materials required:

- Thromboplastin Reagent, e.g. TEClot PT-S (REF A0230)

Storage & Stability

Unopened reagents are stable until the expiration date shown on the label stored at 2°-8°C. Reconstituted reagent:

	2-8 °C	20-25 °C	37°C
Factor II Deficient Plasma	8 hours	4 hours	1 hours

Precautions

Avoid contact with skin and eyes. Wear suitable protective clothing. Dispose components in compliance with local regulations for infectious material. All components are checked for HIV, HBV, HCV. However products from human blood should be considered as potentially infectious.

Specimen collection and storage¹

- Obtain venous blood by clean vein puncture.
- Immediately mix 9 parts blood with 1 part 3.2% sodium citrate (0.105M) and mix well
- Centrifuge the specimen at 1500g for 10 min. (platelet < 10000/µL)
- Separate plasma after centrifugation and store in plastic or siliconised glass tube.
- Use plasma within 4 hours, otherwise store frozen and thaw just prior to use.

Stability of plasma: 4h at 18-26°C 8h at 2-8° 14d at -20°C 6m at -70°C

Procedure

Notes on Procedure:

- Perform all testing in duplicate.
- Prepare dilutions immediately prior to testing.
- A new standard curve must be prepared each time testing is performed.

A. Automated Method

Factor II	A4	A6		A4	A6		A4	A6
PAT Patient	10µL	CP1	8µL	CP1	Incubation	120s	SENS	2
BUF IBS Buffer	40µL	P39	32µL	P79	Maxtime	180s	POINTS	6
CLR -	0µL	-	0µL	-	Unit	32771	MIX	No
DP FII def. Plasma	50µL	P33	40µL	P58	Method	Coag	Clean	0 2
R0 -	0µL	P00	0µL	P00	Math	log XY	Multi	0 2
R1 -	0µL	P00	0µL	P00	CT-Mech	No	S-Corr	0%
R2 PT Reagent	100µL	P25	80µL	P46	Deadtime	9s	T-Corr	0%

B. Manual Method: Coatrom M system

- Pipette **25 µL diluted plasma** (1:10 in IBS buffer) into a test cuvette.
 - Add **25 µL Factor II deficient plasma**
 - Incubate at 37°C for 2 minutes.
 - Add **50 µL Thromboplastin reagent** and simultaneously start test.
- For other instrument, please refer to your instrument manual for more detailed instrument specific instructions.

Calibration

TECal Normal or other commercially reference should be used as calibrator. Prepare and run standards as following:

%Factor activity	Dilution	Plasma	IBS buffer
STD1: 100%	1:10	100µL Plasma	900µL
STD2: 50%	1:20	500µL STD1	500µL
STD3: 25%	1:40	500µL STD2	500µL
STD4: 12.5%	1:80	500µL STD3	500µL
STD5: 6.25%	1:160	500µL STD4	500µL

Plot the clotting times obtained with each standard dilution against % Factor activity using log-log graph paper.

Expected Results

Factor concentration is expressed in IU/mL or % activity (1.00 IU/mL = 100%). Typical normal range is 0.70 – 1.25 IU/mL (70-125%)². However results are influenced by the method of clot detection and can vary from laboratory to laboratory. Each laboratory is recommended to establish its own normal range on the specific instrument used.

Quality Control

TEControl or other commercially control plasma should be used for reliable quality control of performance at a frequency in accordance with good laboratory practice (GLP).

Limitations

Great care must be taken to minimize variations which may occur by seemingly insignificant factors.

A. Specimen Collection. AVOID:

- Use only plastic tubes or siliconised glass.
- Delayed mixing of blood with anticoagulant.
- Contamination with tissue thromboplastin.
- Improper ratio of anticoagulant with blood.
- Hemolyzed, icteric or lipemic samples may interfere optical systems

B. Laboratory Techniques

- Perform tests at 37°C.
- Use only high purity water.
- Optimum pH is 7.0-7.5.

Performance Characteristics

Precision: CV% (within run) CV% (inter-runs)

Normal control < 5.0 < 5.0

Abnormal control < 5.0 < 5.0

(Typical performance on instrument Coatrom M4)

Warranty

This product is warranted to perform in accordance with its labelling and literature. TECO disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event will TECO be liable for any consequential damages arising out of aforesaid express warranty.

References

- National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.
- Lothar Thomas, Labor und Diagnose, 6. Auflage, 2005.

Symbols key:

Expiry date	In Vitro Diagnostic	Biological hazard	Catalogue Number	Consult accompanying documents
Store at 2-8°C	EU conformity	Manufacturer	Lot. Number	Authorized Representative

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B. Manuelle Methode: Coatron M system

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