IVD

REF

P5401-010

Intended Use

This product is intended for the quantitative determination of factor IX in human plasma by clotting assay. Factor IX is a clotting factor which when activated to factor IXa results in the activation of factor X in the presence of factor VIIIa, calcium and phospholipid. Samples with factor IX deficiency have prolonged activated partial thromboplastin time (APTT).

Contents & Determinations

Cat.No.	P5401-010
Factor IX Deficient Plasma*	10x1 mL

Determinations

Coatrom M*	Coatrom A4	Coatrom A6
400	200	250

*contains freeze-dried factor IX depleted human plasma

**Micro method (75µL in total)

Preparation

Reconstitute each vial with 1 mL distilled water. Mix gently and allow standing at room temperature for 20 minutes for complete dissolution.

Further materials required:

- aPTT + CaCl₂ reagent, e.g. TEClot PTT-S (REF A0330)

Storage & Stability

Unopened reagents are stable until the expiration date shown on the label stored at 2°-8°C. Reconstituted reagent:

	2-8 °C	20-25 °C	37 °C
Deficient Plasma	8 hours	4 hours	1 hours

Precautions

Avoid contact with skin and eyes. Wear suitable protective clothing. Dispose components in compliance with local regulations for infectious material. All components are checked for HIV, HBV, HCV. However products from human blood should be considered as potentially infectious.

Specimen collection and storage¹

- Obtain venous blood by clean vein puncture.
- Immediately mix 9 parts blood with 1 part 3.2% sodium citrate (0.105M) and mix well
- Centrifuge the specimen at 1500g for 10 min. (platelet < 10000/µL)
- Separate plasma after centrifugation and store in plastic or siliconised glass tube.
- Use plasma within 4 hours, otherwise store frozen and thaw just prior to use.

Stability of plasma: 4h at 18-26°C 8h at 2-8° 14d at -20°C 6m at -70°C

Procedure**A. Automated Method: Coatrom A**

Factor IX	A4	A6		A4	A6		A4	A6
PAT Patient	10µL	CP1	8µL	CP1	Incubation	180s	SENS	1
BUF IBS Buffer	40µL	P39	32µL	P79	Maxtime	180s	POINTS	6
CLR -	0µL	-	0µL	-	Unit	32771	MIX	No
DP FIX def. Plasma	50µL	P34	40µL	P66	Method	Coag	Clean	0 2
R0 -	0µL	P00	0µL	P00	Math	log XY	Multi	0 2
R1 aPTT reagent	50µL	P31	40µL	P60	CT-Mech	No	S-Corr	0%
R2 CaCl ₂	50µL	P26	40µL	P47	Deadtime	15s	T-Corr	0%

B. Manual Method: Coatrom M

- Pipette **25 µL diluted plasma** (1:10 in IBS buffer) into a test cuvette.
- Add **25 µL Factor IX deficient plasma**
- Incubate at 37°C for 2 minutes.
- Add **25 µL APTT reagent**
- Incubate at 37°C for 3-5 minutes (refers to APTT Reagent kit insert).
- Add **25 µL CaCl₂ reagent** (prewarmed at 37°C) and simultaneously start test.

Calibration

TECal Normal or other commercially reference should be used as calibrator. Prepare and run standards as following:

%Factor activity	Dilution	Plasma	IBS buffer
STD1: 100%	1:10	100µL Plasma	900µL
STD2: 50%	1:20	500µL STD1	500µL
STD3: 25%	1:40	500µL STD2	500µL
STD4: 12.5%	1:80	500µL STD3	500µL
STD5: 6.25%	1:160	500µL STD4	500µL

Plot the clotting times obtained with each standard dilution against % Factor activity using log-log graph paper.

Expected Results

Factor concentration is expressed in IU/mL or % activity (1.00 IU/mL = 100%). Typical normal range is 0.70 – 1.42 IU/mL (70-142%)². However results are influenced by the method of clot detection and can vary from laboratory to laboratory. Each laboratory is recommended to establish its own normal range on the specific instrument used.

Quality Control

TEControl or other commercially control plasma should be used for reliable quality control of performance at a frequency in accordance with good laboratory practice (GLP).

Limitations

Great care must be taken to minimize variations which may occur by seemingly insignificant factors.

A. Specimen Collection, AVOID:

- Use only plastic tubes or siliconised glass.
- Delayed mixing of blood with anticoagulant.
- Contamination with tissue thromboplastin.
- Improper ratio of anticoagulant with blood.
- Hemolyzed, icteric or lipemic samples may interfere optical systems

B. Laboratory Techniques

- Perform tests at 37°C.
- Use only high purity water.
- Optimum pH is 7.0-7.5.

Performance Characteristics

Precision: CV% (within run) CV% (inter-runs)

Normal control < 5.0 < 5.0

Abnormal control < 5.0 < 5.0

(Typical performance on instrument Coatrom M4)

Warranty

This product is warranted to perform in accordance with its labelling and literature. TECO disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event will TECO be liable for any consequential damages arising out of aforesaid express warranty.

References

- National Committee for the National Laboratory (NCCLS) Standards: Collection transport and preparation of blood specimens for coagulation testing and performance of coagulation assays. Document H21-A2, vol. 11, No. 23, 1991.
- Lothar Thomas, Labor und Diagnose, 6. Auflage, 2005.

Symbols key:

	Expiry date		In Vitro Diagnostica		Biological hazard		Catalogue Number		Consult accompanying documents
	Store at 2-8°C		EU conformity		Manufacturer		Lot. Number		Authorized Representative

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	A4	A6
Incubation	180s	
Maxtime	180s	
Unit	32771	
Method	Coag	
Math	log XY	
CT-Mech	No	
Deadtime	15s	

	A4	A6
SENS	1	
POINTS	6	
MIX	No	
Clean	0	2
Multi	0	2
S-Corr	0%	
T-Corr	0%	

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