

Spare part price list for Coatron A6 Plus

valid for year: 2018/2019

TECO

REF	Product name	EUR	USD
18 040 00	Download cable	34,00 €	\$40,80
19 030 01	Reagent adapter (di=22,6)	10,48 €	\$12,60
19 030 03	Reagent adapter (di=24,2)	10,48 €	\$12,60
21 006 01*	Main board Coatron 5000, tested	1.015,00 €	\$1.218,00
21 009 02*	Adapter board	157,50 €	\$189,00
21 010 00	Keyboard	88,00 €	\$105,60
21 011 06	Label, system	24,00 €	\$28,80
21 012 06	Label, reagent	28,00 €	\$33,60
21 021 00	Spring	11,00 €	\$13,20
21 025 00	Heating foil, optic block	36,00 €	\$43,20
21 031 00	Heating foil, reagent block	28,00 €	\$33,60
21 034 00	Waste complete, standard	30,80 €	\$37,00
21 034 50	Waste complete, standard, new series	44,40 €	\$53,28
21 037 00	Emergency fixation	3,20 €	\$3,84
21 048 01	Extender 2, modified (pair)	10,00 €	\$12,00
21 101 02	Laser-Barcodescanner	737,50 €	\$812,00
21 103 00	Standby switch	4,00 €	\$4,80
21 106 02*	Fuse board	210,00 €	\$252,00
21 113 01	Barcodelabel 1-45	9,16 €	\$11,00
21 122 00	Casing left, painted	138,00 €	\$165,60
21 140 00	Extender rack, (pair)	24,00 €	\$28,80
21 157 00	Assay microswitch	18,00 €	\$21,60
21 162 00	Cable optic - data	13,00 €	\$15,60
21 163 00	Cable temp. Sensor	8,00 €	\$9,60
21 163 06	Cable temp. Sensor Reagenzblock 2	8,00 €	\$9,60
21 166 00	Cable reedswitch	10,00 €	\$12,00
21 169 00	Data cable Display	9,64 €	\$11,60
21 170 00	Cable keyboard extension	14,00 €	\$16,80
21 227 00	Reagent Motor cpl.	104,00 €	\$124,80
21 229 00	Shielding display	2,60 €	\$3,12
21 817 06	Optic block cpl. 405 nm, with boards	800,00 €	\$960,00
21 913 00	Assay display cpl.	141,00 €	\$169,20
21 916 00	Assay peltierelement	60,00 €	\$72,00
21 916 06	Assay peltierelement 02	60,00 €	\$72,00
21 921 01	Assay gearboxmotor	150,00 €	\$180,00
21 930 00	Assay sample rack 1	217,50 €	\$261,00
21 931 00	Assay sample rack 2	217,50 €	\$261,00
21 932 00	Assay sample rack 3	217,50 €	\$261,00
23 001 01	Casing, bottom	210,00 €	\$252,00
23 001 05	Casing, rear	300,00 €	\$360,00
23 016 00	Cuvette wheel	165,00 €	\$198,00
23 037 00	Waste complete, for pump kit	64,00 €	\$76,80
23 063 01	Tubing blue (pump - container)	72,00 €	\$86,40
23 063 02	Tubing red	72,00 €	\$86,40
23 063 03	Tubing yellow	72,00 €	\$86,40

23 067 05	Pipette needle	210,00 €	\$252,00
23 070 00	Printer, easy paper load	320,10 €	\$384,12
23 123 00	Drawer cover left, painted	120,00 €	\$144,00
23 124 00	Drawer cover right, painted	135,00 €	\$162,00
23 801 00	Z-Rack with guide	484,00 €	\$580,80
23 803 00	Cable, liquid detection coax L=110	132,00 €	\$158,40
23 804 00	Adapter, tubing mesh to Z rack	60,00 €	\$72,00
23 804 01	Adapter, tubing mesh to Z rack angled (dual Z)	60,00 €	\$72,00
23 805 00	Adapter, Z box to tuding mesh	56,00 €	\$67,20
23 806 00	Support mesh, Z-axis tubing and coax	84,00 €	\$100,80
23 807 00	Home sensor, X-axis	162,00 €	\$194,40
23 808 00	Cover, dual Z-axis	243,00 €	\$291,60
23 810 00	Main board pipetting system	812,50 €	\$975,00
23 811 00	Main board X, Y, Z-axis	812,50 €	\$975,00
23 814 00	PCBA board, dual Z box	625,00 €	\$750,00
23 816 00	Valve XC, teflon M6 with label	390,00 €	\$468,00
23 817 00	Optic block cpl. 405 nm, Coatron 5000	800,00 €	\$960,00
23 818 00	Insulationblock cpl. A6 Plus	446,10 €	\$535,00
23 830 00	Pump with valve complete (23 061 06 + 23 061 08)	2.100,00 €	\$2.520,00