

Intended Use

The Double Cuvette as general laboratory article is suitable to hold a reaction mixture for use in optical analyzers which are used in laboratories for in vitro diagnostic tests.

The cuvettes are ready for immediate use. They have unlimited shelf life if stored at +5-50°C.

Contents

Product	Double Cuvettes
Cat.No.	19 000 02
Content	250 pcs (500 det.)

Cuvettes can be used with Coagulometer Coatron M2/ M4.

Precautions and waste information

The Cuvettes should only be used once in analyzers. To prevent contamination (sample/reagent), it is advised to avoid contact with skin and eyes. Suitable protective clothing and gloves are recommended. Please also note the disposal of components in accordance with local regulations for infectious material.

Warranty

This product is warranted to perform in accordance with its labelling and literature. TECO disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event will TECO be liable for any consequential damages arising out of aforesaid express warranty.



Material: pure, clear Polystyrol

Maximum volume should be less than ~ 500µL (single)

Minimum volume: 75 µL per cuvette (single)

Dimensions max.: (mm) L30 x W11 x H24 (double)
(mm) Ø11 x 19 (single)

Example Picture of the package – 250 Double Cuvettes



Packaging:

1. Log Bag, Dim.: (mm) 180 x 250 x 0,05, Mat.: LDPE;
2. Card Box, Dim.: (mm) 225 x 165 x 60